HOME AND COMMUNITY CARE SUPPORT SERVICES North East

Regional Medical Supplies List Revised July 2022



Table of Contents

| Dressings, Bandages, Tape, Compression & Wraps | 1 |
|---|----------|
| Gauze, Sponges & Packing | 1 |
| Tapes & Securement | 2 |
| Bandages & Wraps Compression Wraps | 3 4 |
| Other, Specialty & Kits | 5 |
| Wound Care | 6 |
| Absorptives | 6 |
| Gelling Fiber & Antimicrobials Silver Gelling Fiber | 6 |
| Alginate & Antimicrobial Silver Alginate Dressings | 7 |
| Antimicrobial & Silver Dressings | 7 9 |
| Composite & Multifunction Dressings Clear Acrylic Dressings | 9 10 |
| Foams | 10 |
| Hydrocolloids | 11 |
| Hydrogel | 11 |
| Films & Membranes Odour Absorbent Dressings | 11 12 |
| Non-adherent | 12 |
| Offloading Devices | 13 |
| Other, Specialty & Kits | 14 |
| Negative Pressure Wound Therapy (NPWT) | 17 |
| Disposable Single-use NPWT | 17 |
| Reusable Foam NPWT | 18 |
| Cleansers & Antiseptic Solutions | 21 |
| Infusion I.V. Supplies | 22 |
| Needles & Syringes | 22 |
| Prefilled Syringes | 23 |
| I.V. Solutions | 23 |
| I.V. Catheters and Access | 24 24 |
| I.V. Gravity Sets & Extension Sets I.V. Securement | 24 25 |
| Infusion Pump & Supplies | 25 25 |
| Blood Collection Supplies | 27 |
| Other IV Supplies & Accessories | 28 |
| I.V. Kits | 29 |
| Ostomy Supplies and Accessories | 38 |
| Ostomy Skin Barriers for 2 Piece System | 38 |
| Urostomy Pouches for 2 Piece System | 39 |
| Drainable Pouches for 2 Piece System | 39 |
| One Piece Urostomy System | 40 |
| One Piece System with Drainable Pouches Pediatric 2 Piece System | 40 41 |
| Ostomy Accessories | 41 |
| Urinary Supplies | 44 |
| Foley Catheters | 44 |
| Intermittent Catheters | 44 |
| Urinary Drainage Bags | 45 |
| Catheter Stabilization | 46 |
| Catheter Accessories | 46 |
| Continuous Bladder Irrigation Urinary Kits | 46 47 |
| Feeding, Suction & Tracheostomy Supplies | 49 |
| Enteral Feeding | 49 49 |
| Suction Catheters, Tubing and Accessories | 49 51 |
| Tracheostomy Products | 51 |
| Oxygen Monitoring | 52 |
| Protective Wear | 53 |

| Appendix A - Depot/Drop-off Locations | А |
|---------------------------------------|----|
| Algoma | A1 |
| Nipissing - Temiskaming | A2 |
| Cochrane | A4 |
| Sudbury-Manitoulin-Parry Sound | A5 |

- The Regional Medical Supplies List identifies the items funded by Home and Community Care Support Services North East that are available to eligible patients who are authorized for Home and Community Care nursing services as per applicable policies.
- If different products are requested, an automatic substitution policy may be initiated to provide the equivalent product as identified within this list.
- The Home and Community Care Support Services North East is not obligated to fund or provide any items that are not identified on this list.
- Home and Community Care Support Services North East-funded nurses are expected to be trained on the supplies identified on this list, and use the supplies in alignment with the manufacturers' recommendations.

Searching the Medical Supplies List Electronically

BOOKMARKS: Make sure the Navigation Pane is showing on the left hand side. If it is not press <F4>

and then click on the second Navigation Pane option that looks like a ribbon. This will display the bookmarks. Click on the bookmark you want to view.

PRODUCTS: Click on the **Edit menu**> **Find**> or press <**Ctrl** + **F** > and type in any word pertaining to product such as product code, description, manufacturer or manufacturer code and then press Enter.

| Find | | × |
|------|----------|------|
| Find | | • |
| | Previous | Next |

To navigate within the document, click the arrows or enter in a specific page number



Please report any errors or feedback to <u>NE.Contracts@hccontario.ca</u>

DRESSINGS, BANDAGES, TAPE, COMPRESSION & WRAPS

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | | |
|---------------------|--|---------------------|----------------------------|-----------------------------------|--------|------------------|---|--|--|
| materials in t | | cut and placed in | a wound because lo | oose fibers (lint) n | nay ge | t lost in the w | to the fibers, wicked away, or passed through into other absorbent yound and delay healing. Nonwoven gauze consists of fibers pressed rbent capacity. | | |
| ROLL BA | NDAGES: Knitted stretch gauze in roll configur | ation used as seco | ondary dressing for | limbs and to hold | dressi | ngs in place. | Includes sterile and non-sterile. | | |
| SDR-0020 | Sterile Small Woven Gauze 2" Roll | Gazetex | Derma Sciences DUP9321 | 2.25" x 3 yards 5.7 cm x 2.7 m | Y | 7 | Gazetex bandage rolls are washed, low-linting woven gauze rolls; 6- ply. These bandages are used for wrapping or packing large and deep wounds. When used as a primary dressing, the dressing provides fast | | |
| SDR-0021 | Sterile Large Woven Gauze 4" Roll | Gazetex | Derma Sciences DUP9322 | 4.5" x 4 yards 11.4cm x 3.7m | Y | 7 | wicking, absorbency, and aeration. When used as a secondary dressing, the dressing provides additional cushioning and comfort. Does not contain natural rubber latex | | |
| SDR-0007 | Conforming Stretch Bandage Wrap 2" (Roll) | Duform | Derma Sciences DUP75102 | 2" x 4.1 yards 5 cm x 3.75 m | Y | 7 | Synthetic knitted stretch bandage ideal as a secondary dressing. N | | |
| SDR-0008 | Conforming Stretch Bandage Wrap 4" (Roll) | Duform | Derma Sciences DUP75104 | 4"x 4.1 yards 10cm x 3.75m | Y | 7 | Sterile | | |
| SDR-0060 | Conforming Cotton Bandage 4" (Roll) | Dutex | Derma Sciences DUP76783 | 4" x 4.1 yards 10cm x 3.75m | Y | 7 | Cotton, 2 ply woven bandage ideal as a secondary dressing. Non- stretch and Non-sterile | | |
| DRESSIN | GS AND PACKING: Sterile dressing product | s to cover surgical | incisions to control | bleeding, absorb | exuda | ite, provide p | rotection, and non-sterile products to secure primary dressings in place. | | |
| SDR-0013 | Oval Eye Pad Non-woven Sterile | Surecare | Derma Sciences DUP82911 | 1-5/8" x 2-5/8" 4 x 6.5 cm | Y | 7 | | | |
| SDR-0001 | Abdominal Pad Sterile 5"x9" | Dupad | Derma Sciences DUP87059 | 5" x 9" 12.7 x 22.8 cm | Y | 7 | Non-woven; non-adherent | | |
| SDR-0002 | Abdominal Pad Sterile 8"x10" | Dupad | Derma Sciences DUP87810 | 8" x 10" 20 x 25 cm | Y | 7 | Non-woven; non-adherent | | |
| SDR-0003 | Gauze Sterile Woven 2"x2" (2 per Pack) | Ducare | Derma Sciences DUP82212 | 2" x 2" 5 x 5 cm | Y | 14 Packs of 2 | Cotton, 12 Ply, Pack of 2 | | |
| SDR-0004 | Gauze Sterile Woven 4"x4" (2 per Pack) | Ducare | Derma Sciences DUP82408 | 4" x 4" 10 x 10cm | Y | 14 Packs of 2 | Cotton, 8 Ply, Pack of 2 | | |
| SDR-0017 | Sterile Packing Plain 1/4" (Bottle) | PAK-ITS | Derma Sciences DUP59120 | 1/4" x 5 yards 0.64cm x 4.5m | Y | 2 | Single-ply ribbon gauze strips that can be used dry to absorb drainage | | |
| SDR-0018 | Sterile Packing Plain 1/2" (Bottle) | PAK-ITS | Derma Sciences DUP59220 | 1/2" x 5 yards 1.27cm x 4.5m | Y | 2 | in exuding wounds or soaked with a medication of choice. A secondary dressing is needed to hold the strips in place and to protect | | |
| SDR-0019 | Sterile Packing Plain 1" (Bottle) | PAK-ITS | Derma Sciences DUP59320 | 1" x 5 yards 2.54cm x 4.5m | Y | 2 | the wound from outside contaminants. (Bottles not for use at clinic - See single use packs below). | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes |
|---------------------|--|---------------|-----------------------------|---------------------------|---------|------------------|--|
| SDR-0022 | Sterile Gauze Packing Cotton 2x36" | PAK-ITS | Derma Sciences DUP564036 | 2" x 36" 5 cm x 90 cm | Υ | 3 | |
| SDR-0088 | Sterile Gauze Packing Cotton Single Use 1/4" | PAK-ITS | Derma Sciences DUP561036 | 1/4" x 36" 0.64 x 90cm | Υ | 3 | |
| SDR-0089 | Sterile Gauze Packing Cotton Single Use 1/2" | PAK-ITS | Derma Sciences DUP4853 | 1/2" x 54" 1.3 x 137cm | Y | 3 | |
| SDR-0090 | Sterile Gauze Packing Cotton Single Use 2" | PAK-ITS | Derma Sciences DUP59420 | 2" x 54" 5 cm x 4.5m | Υ | 3 | |
| SDR-0023 | Sterile Gauze Packing Cotton Fan Folded 2x72" | PAK-ITS | Derma Sciences DUP4869 | 2" x 72" 5 cm x 182 cm | Y | 3 | |
| SDR-0064 | Раск) | DuSoft | Derma Sciences DUP84122 | 2" x 2" 5 x 5 cm | Υ | 14 Packs of 2 | Synthetic, 4 Ply, Pack of 2 |
| SDR-0065 | Sponge Synthetic Sterile Non-Woven 4"x4" (2 per Pack) | DuSoft | Derma Sciences DUP84144 | 4" x 4" 10 x 10cm | Υ | 14 Packs of 2 | Synthetic, 4 Ply, Pack of 2 |
| SDR-0076 | Sponge Synthetic Non-sterile Non-Woven 4"x4" | DuSoft | Derma Sciences DUP94144 | 4" x 4" 10.2 x 10.2 cm | Υ | 1 Pack of 200 | Synthetic, 4 Ply, Pack of 200 |
| Available in a | | The non-woven | sponges provide su | • | | | ng around a tracheostomy tube or any type of tube, IV line/PICC line. fast-wicking, and non-adherent to wound beds. The woven gauze |
| | Drain Gauze Sponge Tracheostomy Sterile Woven 4"x4" | Ducare | Derma Sciences DUP82910 | 4" x 4" 10 x 10cm | Y | 14 | Woven, 12 Ply, Pack of 1 |
| SDR-0058 | Drain Sponge Sterile Non-woven 4x4 | DuSoft | Derma Sciences DUP84916 | 4" x 4" 10 x 10cm | Y | 14 Packs of 2 | 6 Ply Non-woven, Pack of 2 |
| SOT-0036 | Drain Sponge Sterile Non-woven 2x2 | DuSoft | Derma Sciences DUP84918 | 2" x 2" 5 x 5 cm | Υ | 14 Packs of 2 | 6 Ply Non-woven, Pack of 2 |
| SOT-0024 | Drain Sponge Antimicrobial Sterile Non-woven 4x4 | Excilon | Covidien Z7088 | 4" x 4" 10 x 10cm | Y | 14 Packs of 2 | PHMB, 6 Ply Non-woven, Pack of 2 |
| | SECUREMENT: Medical and surgical tapes us | | ent aid in dressing o | or tubings (NG, I\ | /, Fole | y, etc.). | |
| SDR-0039 | Non Allergenic Self-Adhesive Dressing Retention Tape 2" | Ultrafix | Derma Sciences DUP62036 | 2" x 1 yrd 5 cm x 1 m | Y | /m | A non-sterile and non-woven fixative tape that is breathable and gentle to the skin. Custom cut by metre lengths |
| SDR-0040 | Non Allergenic Self-Adhesive Dressing Retention Tape 4" | Ultrafix | Derma Sciences DUP62037 | 4" x 1 yrd 10 cm x 1 m | Υ | / m | A non-sterile and non-woven fixative tape that is breathable and gentle to the skin. Custom cut by metre lengths |
| SDR-0048 | Montgomery Straps and Ties | Montgomery | Derma Sciences DUP97991 | 18 cm x 28 cm | Y | 1 | Non-sterile, non-woven |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
|---------------------|--|--------------------------------------|-------------------------------|---------------------------------|-----|---------------|---|--|
| SDR-0086 | White Surgical Tape 1" Single Use Roll | Transpore (or equivalent) | 3M 1534S-1 (or equivalent) | 1" x 1.5 yard 2.5 cm x 1.3 m | Υ | 1 | Universal hypoallergenic tape can be used for all applications. | |
| SDR-0087 | Blue Silicone Tape 1" Single Use Roll | Micropore S | 3M 2770S-1 (or equivalent) | 1" x 54" 2.5 cm x 1.3 m | Υ | 1 | Only order for patients with very fragile at risk skin who react to SDR-0086 1" White Surgical Tape | |
| SDR-0037 | Wound Closure Strip 1/2" (6 Strips) Sterile | Primaflex | priMED Medical PM3-71241 | 1/2" x 4" 12 x 100 mm | Υ | 2 Packs | 6 Strips per pack | |
| SDR-0038 | Wound Closure Strip 1/4" (3 Strips) Sterile | Primaflex | priMED Medical PM3-71431 | 1/4" x 3" 6 x 75 mm | Υ | 2 Packs | 3 Strips per pack | |
| cleansing, pa | cking or to provide absorption as a second layer or clude securing IVs and splints and providing mild o | dressing, to hole compression/sup | d a wound cover in | place, to apply pr | | | sition may include cotton, polyester, rayon or nylon. They function for n a wound. Some are sterile or non-sterile in bulk and may be adhesive. | |
| SDR-0027 | Tubular Net Elastic Bandage - Fingers/Toes/Small Wrist | Medline Elastic Net | NONNET02 | 7 2/3" 19.5 cm | Y | 2 x 1 metre | | |
| SDR-0083 | Tubular Net Elastic Bandage - Wrist & Elbow | Medline Elastic Net | NONNET03 | 9 1/3" 23.7 cm | Y | 2 x 1 metre | | |
| SDR-0028 | Tubular Net Elastic Bandage - Elbow & Arm | Medline Elastic Net | NONNET04 | 11 2/3" 29.6 cm | Y | 2 x 1 metre | | |
| SDR-0029 | Tubular Net Elastic Bandage - Knee & Thigh | Medline Elastic Net | NONNET05 | 12 3/8" 31.4 cm | Y | 2 x 1 metre | Non-adhesive and holds dressing/devices securely in place while allowing freedom of movement. Roll of 25 yards cut to 1 meter lengths | |
| SDR-0084 | Tubular Net Elastic Bandage - Shoulder Hip Thigh | Medline Elastic Net | NONNET06 | 17" 61 cm | Y | 2 x 1 metre | | |
| SDR-0030 | Tubular Net Elastic Bandage - Chest Back Perineum Axilla | Medline Elastic Net | NONNET07 | 28 2/3" 72.8 cm | Y | 2 x 1 metre | | |
| SDR-0085 | Tubular Net Elastic Bandage - Chest Region Large | Medline Elastic Net | NONNET08 | 31 1/3" 79.6 cm | Y | 2 x 1 metre | | |
| SDR-0031 | Tubular Net Elastic Bandage - Chest Region XL | Medline Elastic Net | NONNET09 | 35" 88.9 cm | Y | 2 x 1 metre | Non-adhesive and holds dressing/devices securely in place while allowing freedom of movement. SPECIAL ORDER ITEM - May take u to 5 days for stock. | |
| SDR-0035 | Tubular Net Elastic Bandage - Chest Region XXL | Medline Elastic Net | NONNET10 | 41" 93.1 cm | Y | 2 x 1 metre | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
|---------------------|---|------------------|------------------------------|----------------------|---------|---------------|---|--|
| SCM-0003 | Mesh Pants XL / XXL | Wings | Cardinal Health Z707A | 40"-56" | Υ | 1 | Incontinence supplies are not funded as per policy. These mesh pants | |
| SCM-0004 | Mesh Pants Large / XL | Wings | Cardinal Health Z706A | 32"-42" | Y | 1 | are only authorized for patients requiring this product as a means of containing wound care dressings to coccyx, perineal area, etc. | |
| SCM-0005 | Mesh Pants Small / Medium | Wings | Cardinal Health Z705A | 22"-32" | Y | 1 | oomanning wound ouro aroooningo to ooooyx, pormour aroa, oto. | |
| SDR-0068 | Multipurpose Tubular Elastic Bandage Size C | Reliamed | Cardinal Health ZG275TB | 2.75" x 36" | Y | 2 | - Limb sizes C through G available. Other sizes may available only as | |
| SDR-0069 | Multipurpose Tubular Elastic Bandage Size D | Reliamed | Cardinal Health ZG30TB | 3" x 36" | Y | 2 | a exception item special order with approval. Sizes intended as a guideline only and should be used in conjunction with existing clinical | |
| SDR-0070 | Multipurpose Tubular Elastic Bandage Size E | Reliamed | Cardinal Health ZG35TB | 3.5" x 36" | Y | 2 | protocols that consider physiological state and wound and skin condition. | |
| SDR-0071 | Multipurpose Tubular Elastic Bandage Size F | Reliamed | Cardinal Health ZG40TB | 4" x 36" | Y | 2 | - Order in 1 yard lengths (36") - If overlapped, applies pressure of 10-15 mm Hg | |
| SDR-0072 | Multipurpose Tubular Elastic Bandage Size G | Reliamed | Cardinal Health ZG45TB | 4.5" x 36" | Y | 2 | Can be washed (gentle cycle & air dry) and reused multiple times Contains latex | |
| SVS-0003 | High Stretch Elastic Bandage Non-Latex 4" | Dusor | Derma Sciences DUP71414LF | 4" x 5 yards | Υ | 1* | *May order 2 for initial order only if needed. Non-Sterile | |
| SVS-0004 | High Stretch Elastic Bandage Non Latex 6" | Dusor | Derma Sciences DUP71416LF | 6" x 5 yards | Y | 1* | *May order 2 for initial order only if needed. Non-Sterile | |
| COMPRES | SION WRAPS: Short Stretch multi-layer elast | c stretch bandag | e system with limite | ed extensibility (le | ss than | າ 100%) for ເ | se to provide moderate to high compression | |
| SVS-0006 | 2 Layer Compression System 35-40mmHg | Coban 2 | 3M 2094 | Each | Y | 2 | 2 Rolls - Comfort layer and compression layer. Ideal for majority of patients with venous leg ulcers, lymphedema and other conditions where compression bandage therapy is appropriate. | |
| SVS-0007 | 2 Layer Compression System 25-30 mmHg | Coban 2 Lite | 3M 2794 | Each | Y | 2 | 2 Rolls - Comfort layer and compression layer. Designed to be more comfortable for patients less tolerant of compression bandages. | |
| SVS-0009 | Tubular Elastic Compression Small Navy Stripe | EdemaWear | CA600001 | Each Pair (2) | Y | 1/Month | Up to 2 Pair at admission then 1 Pair per month as needed | |
| SVS-0010 | Tubular Elastic Compression Medium Yellow Stripe | EdemaWear | CA960001 | Each (1) | Y | 1/Month | 15-20 mmHg; Air-permeable; Can be cut Small (Navy): Up to 18″/45cm circumference (Foot/Calf/Arm) Medium (Yellow): Up to 30″/75cm circumference (Knee/Thigh/Arm) | |
| SVS-0011 | Tubular Elastic Compression Large Red Stripe | EdemaWear | B120L01 | Each (1) | Y | 1/Month | Large (Red): Up to 46"/115cm circumference (Thigh/Torso) | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | | |
|---------------------|---|-----------------|-----------------------|-----------|-----|---------------|---|--|--|
| OTHER & | SPECIALTY: | | | | | | | | |
| SCM-0019 | Blue Absorbent Pad 23"x24" | Cardinal Health | Cardinal UPSMD2324 | 23" x 24" | Y | 1/Visit | Incontinance supplies are not funded as per approved policy. This product is authorized only for patients for containment of drainage or for methotrexate administration in-home/clinic ONLY. Approved clinics may order in bulk. | | |
| DRESSIN | G KITS: | | | | | | | | |
| SDR-0005 | Dressing Tray Kit - Disposable Single Use | MED-Rx | CHS 85-4009 | Each Kit | Y | 1 | For catheter treatments only (PNT or PICC). Not for use in wound care. | | |
| 3 Compa | artment 1200 CC Outer Tray | | | 1 | | | | | |
| Rayon C | otton Balls Large Absorbent | | | 6 | | | | | |
| Gauze 4' | " x 4" - 8 Ply | | | 4 | | | | | |
| 4" Plastic | c Transfer Forcep | | | 3 | | | | | |
| Absorber | Absorbent Towel - 3 Ply | | | | | | | | |
| CSR Ste | rile Field Towel | | | 1 | | | | | |
| Waste Ba | ag with Tie | | | 1 | | | | | |

WOUND CARE

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
|---------------------|-----------------------------------|----------------------------|----------------------|------------------------|-----|---------------|---|-----------------------|
| Home and | Communuty Care Support Services I | - | | - | | | on all wound care supplies that are funded. No for use, contraindications, wear time, etc. | irses are expected to |
| | | | | | | | phly absorptive layers of fibers, such as cellulose, cotton incisions, lacerations, abrasions, burns, donor or skin gra | |
| SWC-0107 | Absorbent Dressing 4"x6" | Mesorb | Mölnlycke 677201 | 4" x 6" 10 x 15 cm | Y | 7 | | 1 |
| SWC-0108 | Absorbent Dressing 6"x8" | Mesorb | Mölnlycke 677601 | 6" x 8" 15 x 20 cm | Y | 7 | | |
| SWC-0109 | Absorbent Dressing 8"x12" | Mesorb | Mölnlycke 677801 | 8" x 12" 20 x 30 cm | Y | 7 | | |
| SWC-0142 | Super-Absorbent Dressing 4"x6" | Mextra Super- absorbent | Mölnlycke 610100 | 4" x 6" 10 x 15 cm | Y | 4 | - For best result, the pad should overlap the dry | |
| SWC-0143 | Super-Absorbent Dressing 6"x8" | Mextra Super- absorbent | Mölnlycke 610300 | 6" x 8" 15 x 20 cm | Y | 4 | surrounding skin by at least 2 cm. May be left in place for several days. The dressing should be changed according to the clinical condition | |
| SWC-0144 | Super-Absorbent Dressing 8"x12" | Mextra Super- absorbent | Mölnlycke 610500 | 8" x 12" 20 x 30 cm | Y | 4 | of the wound, or when saturated. | |
| Converts to a | | | | | | | e of action" to absorb wound exudate and maintain a mo operties • Supports autolytic debridement • Best for mod | |

| | - | | | | | | | |
|----------|---|---------------------|--------------------|------------------------|---|---|--------------|---|
| SWC-0150 | Gelling Fiber Dressing 4"x4" | Aquacel Extra | ConvaTec 420672 | 4" x 4" 10 x 10 cm | Y | 4 | | |
| SWC-0151 | Gelling Fiber Dressing 6"x6" | Aquacel Extra | ConvaTec 420673 | 6" x 6" 15 x 15 cm | Y | 4 | | |
| SWC-0106 | Gelling Fiber Ribbon 1"x18" | Aquacel | ConvaTec 403770 | 1" x 18" 2 x 45 cm | Y | 4 | | |
| SWC-0152 | Gelling Fiber Antimicrobial Silver Dressing 4"x4" | Aquacel Ag Extra | ConvaTec 420676 | 4" x 4" 10 x 10 cm | Y | 4 | | 8 |
| SWC-0153 | Gelling Fiber Antimicrobial Silver Dressing 8"x12" | Aquacel Ag Extra | ConvaTec 420679 | 8" x 12" 19 x 29 cm | Y | 4 | | |
| SWC-0032 | Gelling Fiber Antimicrobial Silver Impregnated Ribbon 1"x18" | Aquacel Ag | ConvaTec 403771 | 1" x 18" 2 x 45 cm | Y | 4 | Ionic Silver | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | ΑΡΟ | Weekly Max | Notes |
|---------------------|---|--------------------|----------------------------|------------------------------|-----|---------------|--|
| exudate, suc | | ers and venous ir | nsufficiency ulcers • | Has hemostatic | | | nposed of seaweed derivatives • Indicated for wounds with moderate to heavy quires a secondary dressing • Should not be used on dry wounds • Lower tensile |
| SWC-0036 | Alginate Dressing 4"x4" | Tegaderm | 3M 90112 | 4" x 4" 10 x 10 cm | Υ | 3 | |
| SWC-0120 | Antimicrobial Silver Alginate Dressing | Silvercel | KCI Medical CAD011 | 4.3" x 4.3" 11 x 11 cm | Y | 3 | |
| SWC-0154 | Antimicrobial Silver Rope Alginate Dressing | Silvercel | KCI Medical CAD7230 | 2.5 x 30.5cm | Y | 3 | |
| sites and sur | | ty to handle exua | | | | | e risk of infection in partial- and full-thickness wounds, over percutaneous line of the product. Does not replace systemic antibiotics for deeper tissue infections. |
| SWC-0034 | Non-adherent Antiseptic Impregnated Gauze 4"x4" | Bactigras | Smith & Nephew 66003650 | 4" x 4" 10 x 10 cm | Y | 4 | For adjunctive treatment and prevention of infection in skin loss lesions, including wounds, burns and ulcers |
| SSP-0003 | Ointment with lodine 10Gr Tube | lodasorb | Smith & Nephew 66060630 | 10 Gram Tube | Y | 1 | Cadexomer lodine should not be used in patients with known or suspected lodine sensitivity, history of any thyroid disorder, or pregnant or lactating women |
| SDR-0078 | Non-adhesive Povidone-Iodine Dressing 2"x2" | Inadine | KCI Medical P01481 | 2" x 2" 5 x 5 cm | Y | 3 | 10% Povidone lodine; equivalent to 1.0% available iodine |
| SDR-0079 | Non-adhesive Povidone-Iodine Dressing 3.75"x3.75" | Inadine | KCI Medical P01512 | 3.75 x 3.75" 9.5 x 9.5 cm | Y | 3 | 10% Povidone lodine; equivalent to 1.0% available iodine |
| SWC-0119 | Silver Antimicrobial Barrier Dressing Non- Absorbent 4"x4" | Acticoat Flex 3 | Smith & Nephew 66800399 | 4" x 4" 10 x 10cm | Y | 3 | Not compatible with oil-based products, such as petrolatum |
| SWC-0156 | Antimicrobial Silver Mesh | Acticoat Flex 7 | Smith & Nephew 66800397 | 4" x 5" 10 x 12.5cm | Y | 1 | Up to 7 day wear time. Not to be ordered with dressing changes more frequently than Q3 days. Can be cut to fit. |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
|---------------------|---|-------------------------|--|----------------------------------|-----|----------------|---|--|
| SWC-0010 | PHMB Antimicrobial Gauze Square 4"x4" | Curity | Cardinal Health Z2533 | 4" x 4" 10.2 x 10.2 cm | Y | 7 | Impregnated with 0.2% PHMB; 12 Ply | |
| SWC-0112 | PHMB Antimicrobial Non-Woven Sponges 4"x4" | Kendall | Cardinal Health Z8044AMD or equivalent | 4" x 4" 10.1 x 10.1 cm | Y | 7 | Impregnated with 0.2% PHMB | |
| SWC-0113 | PHMB Antimicrobial Gauze Roll 4.5" | Kerlix | Cardinal Health Z3332 | 4.5" x 4 1/8 yd 11.4cm x 3.7m | Y | 7 | Impregnated with 0.2% PHMB; 6 Ply | |
| SWC-0115 | Antimicrobial Gauze PHMB Ribbon 1/4" x 1 yard | Curity | Cardinal Health Z7831AMD | 1/4" x 36" 0.6 x 91.4 cm | Y | 7 | Impregnated with 0.2% PHMB | |
| SWC-0116 | Antimicrobial Gauze PHMB Ribbon 1/2" x 1 yard | Curity | Cardinal Health Z7832AMD | 1/2" x 36" 1.3 x 91.4 cm | Y | 7 | Impregnated with 0.2% PHMB | |
| SWC-0117 | Antimicrobial Gauze PHMB Ribbon 1" x 1 yard | Curity | Cardinal Health Z7833AMD | 1" x 36" 2.5 x 91.4 cm | Y | 7 | Impregnated with 0.2% PHMB | |
| SWC-0118 | Antimicrobial Gentian Violet Methylene Blue Dressing 4x4" | Hydrofera Blue | Hydrofera HB4414 | 4" x 4" 10 x 10 cm | Y | 3 | | |
| SWC-0137 | Antimicrobial Gentian Violet Methylene Blue Dressing 6x6" | Hydrofera Blue | Hydrofera HB6614 | 6" x 6" 15 x 15 cm | Y | 3 | - Not to be used longer than 6 months without consuting Wound Specialty Nurse. | |
| SWC-0155 | Antimicrobial Gentian Violet Methylene Blue Moist Foam Dressing 4x5" | Hydrofera Blue Ready | Hydrofera HBRS4520 | 4" x 5" 10 x 12.5cm | Y | 1 | - Not to be used on pregnant or nursing patients. | |
| SWC-0138 | Antimicrobial Gentian Violet Methylene Blue Tunneling Dressing 6x6" | Hydrofera Blue | Hydrofera HBT0906 | 9 mm (Pack of 2) | Y | 3 | | |
| SWC-0131 | Antimicrobial Gel with Silver 0.25 oz | SilvaSorb | Medline MSC93025EP | 0.25 oz Tube | Y | As required | For clinic patients only | SILVASORB' GEL Bayes Antibilization of the second s |
| SWC-0132 | Antimicrobial Gel with Silver 1.5 oz | SilvaSorb | Medline MSC9301EP | 1.5 oz Tube | Y | 1 | Community patients - not clinic | Andrewing Stranger Andrewing Stranger Andrew |
| SWC-0134 | Textile Antimicrobial Silver (Custom order by 12" lengths) | InterDry AG | Coloplast 7910 | 10" x Custom 12" lengths | Y | 4 | Moisture-wicking fabric. Custom order per 12" lengths | <image/> <section-header></section-header> |
| SWC-0157 | Foam Silver Dressing | Biatain AG | Coloplast 9622 | 4" x 4" 10 x 10 cm | Y | 4 | Non-adhesive | Annen An Annen An Annen An Annen An |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes |
|---------------------|--|-----------------------|----------------------|---------------------------------|-----|---------------|---|
| SWC-0145 | Silver Non-Adhesive Pad Dressing | PolyMem Silver | 1044 | 4.25" X 4.25" 10.8 x 10.8 cm | Y | 2 | Cleanses, fills, absorbs and moistens wounds. Combined primary/secondary dressing or as a secondary dressing. Can be cut to size. |
| SWC-0146 | Silver Non-Adhesive Pad Dressing (MAX) | PolyMem MAX Silver | 1045 | 4" X 4" 10 x 10 cm | Y | | MAX is a thicker PolyMem configuration for wounds with heavier drainage or when longer wear time is desired. |

COMPOSITE & MULTIFUNCION DRESSINGS: Multilayered, combination dressings to provide multiple functions, such as a bacterial barrier, absorption and adhesion. Usually, they are comprised of multiple layers and incorporate a semi- or non-adherent pad that covers the wound. Some are appropriate for autolysis. May also include an adhesive border of non-woven fabric tape or transparent film. Composite dressings can function as either a primary or a secondary dressing on a wide variety of wounds.

| SWC-0147 | Non-Adhesive Pad Dressing | PolyMem | 5044 | 4" x 4" 10 x 10 cm | Y | 2 | Cleanses, fills, absorbs and moistens wounds. Combined primary/secondary dressing or as a secondary dressing. Can be cut to size. | |
|----------|---|---------------------|---------------------------|---------------------------------|---|---|--|--|
| SWC-0148 | Non-Adhesive Pad Dressing (MAX) | PolyMem MAX | 5045 | 4.25" x 4.25" 10.8 x 10.8 cm | Y | 2 | - MAX is a thicker PolyMem configuration for wounds with heavier drainage or when longer wear time is desired. | |
| SWC-0090 | Protease Inhibitor Dressing with Collagen | Promogran | KCI M772028 | 28 cm ² Hexagonal | Y | 2 | Review manufacturer's information regarding contraindications. | |
| SWC-0158 | Protease Inhibitor Dressing with Silver | Promogran Prisma | KCI PS2028 | 28 cm ² Hexagonal | Y | 2 | Protease inhibitor should be discontinued if there is no progression towards wound healing (less than 10% healing in 3 weeks). | |
| SVS-0008 | Non-Elastic Zinc Paste Bandage | Primer | Derma Sciences GL300-1 | 7.6 cm x 9 m | Y | 2 | A gauze bandage that has been evenly impregnated with a non-hardening paste of zinc oxide. Ideal for managing: venous ulcers, lymphaticedema, minor sprains, strains, fractures and thrombophlebitis. Contains no gelatin which prevents the hardening or caking of the bandage. Made without preservatives and with non-raveling thread locked gauze. | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes |
|--------------------------------|---|-----------------------------------|----------------------|-------------------------------------|---------------|---------------|---|
| without chang | | Extended wear | | | | | rproof film. Enables clinicians to monitor small to moderately exuding wounds ial lacerations • Post-operative wounds • First (superficial thickness) and second |
| SWC-0139 | Absorbent Clear Acrylic Small Oval 7.6x9.5cm | Tegaderm | 3M 90800 | 3 in. x 3-3/4 in. 7.6 x 9.5 cm | Y | 1 | For pressure ulcers, skin tears and other low to |
| SWC-0140 | Absorbent Clear Acrylic Medium Oval 11x12cm | Tegaderm | 3M 90801 | 4-3/8 in. x 5 in. 11.1 x 12.7 cm | Y | 1 | - moderately exuding wounds. Dressing should overlap the wound margins by at least 3/8" (1 cm). |
| combination v a transparent | with other materials. Absorption capability depends | s on thickness an | d composition. The | area in contact v | , vith the | e wound s | en cells capable of holding fluids. They may be impregnated or layered in surface is non-adhesive for easy removal. Available with an adhesive border and/or -exudating wounds, depending on product. Occlusive foams without silver should |
| SWC-0099 | Foam Non-Border 4"x4" | Mepilex XT | Mölnlycke 211100 | 4" x 4" 10 x 10 cm | Y | 4 | |
| SWC-0100 | Foam Non-Border 8"x8" | Mepilex XT | Mölnlycke 211400 | 8" x 8" 20 x 20 cm | Y | 4 | |
| SWC-0101 | Foam Non-Border Soft Silicone 6"x6" | Mepilex Lite | Mölnlycke 284300 | 6" x 6" 15 x 15 cm | Y | 4 | |
| SWC-0102 | Foam Border Silicone 3"x3" | Mepilex Border | Mölnlycke 295200 | 3 "x 3" 7.5 x 7.5 cm | Y | 4 | |
| SWC-0103 | Foam Border Silicone 4"x4" | Mepilex Border | Mölnlycke 295300 | 4" x 4" 10 x 10 cm | Y | 4 | |
| SWC-0104 | Foam Border Silicone 6"x6" | Mepilex Border | Mölnlycke 295400 | 6" x 6" 15 x 15 cm | Y | 4 | |
| SWC-0105 | Foam Border Silicone 4"x4" | Mepilex Border Lite | Mölnlycke 281300 | 4" x 4" 10 x 10 cm | Y | 4 | Safeta C |
| SWC-0129 | Foam Non-Adhesive with Ibuprofen 4"x4" | Biatain IBU | Coloplast 4110 | 4" x 4" 10 x 10 cm | Y | 4 | And the second se |
| SWC-0130 | Foam Non-Adhesive with Ibuprofen 6"x6" | Biatain IBU | Coloplast 4115 | 6" x 6" 15 x 15 cm | Y | 4 | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes |
|-------------------------------|---|--|--|------------------------------------|--------------------|---------------|---|
| | | d of gelatin nectir | | olluloso Moistur | o roton | - | ing, contributes to autolytic debridement. Creates occlusive barrier against |
| | · · · · · · | • | | | | | eels and sacral ulcers. Indicated for partial and full thickness wounds with or |
| | | | | | | | hen used on fragile skin and should not be used on heavily draining or infected |
| ounds. | | 5 5 | | | | - | 5 5 5 |
| WC-0051 | Hydrocolloid Regular Square 4x4" | Comfeel Plus | Coloplast 33110 | 4" x 4" 10 x 10 cm | Y | 3 | allower and the |
| WC-0052 | Hydrocolloid Regular Square 6x6" | Comfeel Plus | Coloplast 33115 | 6" x 6" 15 x 15 cm | Y | 3 | Street Street |
| WC-0080 | Hydrocolloid Regular 8"x8" | Comfeel Plus | Coloplast 33120 | 8" x 8" 20 x 20 cm | Y | 3 | Constant Constant Constant |
| WC-0056 | Hydrocolloid Regular Sacral | Comfeel Plus | Coloplast 33285 | 7" x 8" 18 x 20 cm | Y | 3 | |
| WC-0045 | Hydrocolloid Transparent Square 4x4" | Comfeel Plus Transparent | Coloplast 33533 | 4" x 4" 10 x 10 cm | Y | 3 | |
| WC-0048 | Hydrocolloid Transparent Square 6x6" | Comfeel Plus Transparent | Coloplast 33539 | 6" x 6" 15 x 15 cm | Y | 3 | |
| WC-0047 | Hydrocolloid Transparent Rectangle 3.5x5.5" | Comfeel Plus Transparent | Coloplast 33536 | 3.5" x 5.5" 9 cm x 14 cm | Y | 3 | |
| SWC-0053 | Hydrocolloid Transparent Oval | Comfeel Plus Transparent | Coloplast 33283 | 3.5" x 4" 9 x 11 cm | Y | 3 | |
| ehydrate wo | | edded into gauze. | Can decrease pair | | | | and maintain a moist healing environment. The high moisture content serves to ickness wounds, wounds with necrosis, minor burns and radiation tissue dam |
| WC-0023 | Hydrogel (Low Pectin, Low Saline) | Intrasite Gel | Smith & Nephew 7311 | 15 g | Y | 3 | |
| noisture vapo ndicated for | our and atmospheric gases. Moisture vapour tran | nsmission rate (M) wounds with necr | /TR) varies from fill osis, and as both a | m to film. Can he primary or secor | elp red ndary c | uce frictio | Indhesive. They are impermeable to liquid, water and bacteria but permeable to on to susceptible skin (e.g. heels). Transparency allows visualization of the wo Also used to cover IV sites, donor sites, lacerations, abrasions and second-deq |
| | Transparent Film Dressing 10x12cm | Tegaderm | 3M | 4" x 4.75" | Y | 3 | |

| SWC-0072 | Transparent Film Dressing 10x12cm | Tegaderm | 3M 1626W | 4" x 4.75" 10 x 12 cm | Y | 3 | | |
|----------|-------------------------------------|----------|-------------|--------------------------|---|------------|---|--|
| SWC-0128 | Transparent Film Dressing 15x20cm | Tegaderm | 3M 1628 | 6" x 8" 15 x 20 cm | Y | 3 | | |
| SWC-0159 | Transparent Film Dressing 10cm x 1m | Tegaderm | 3M 16004 | 10 cm x 1 m | Y | 1 Meter | Roll cut in 1 meter lengths per patient order | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
|---------------------|--|-----------------------------|--------------------------|------------------------------|--------|---------------|--|----------|
| ODOUR A | BSORBENT DRESSINGS: Odour absorbent | charcoal contair | ned within product t | hat masks the od | our. E | nsure tha | t dressing edges are sealed to control odour and treat the cause if p | ossible. |
| SWC-0037 | Activated Charcoal Dressing with Silver 4"x4" | Actisorb Plus Silver 220 | KCI Systagenix MAS105 | 4" x 4" 10.5 x 10.5 cm | Y | 2 | Silver is not released into the wound, therefore can | |
| SWC-0135 | Activated Charcoal Dressing with Silver 2.5 x 3.75" | Actisorb Plus Silver 220 | KCI Systagenix MAS065 | 2.5" x 3.75" 6.5 x 9.5 cm | Y | 2 | be use on silver-sensitive patients. Remains effective to reduce odour when wet or dry and can be left in place for up to 7 days. NOTE: This product | |
| SWC-0136 | Activated Charcoal Dressing with Silver 4" x 7.5" | Actisorb Plus Silver 220 | KCI Systagenix MAS190 | 4" x 7.5" 10.5 x 19 cm | Y | 2 | should not be cut | |
| | ERENT: Porous sheets of dressings with low adl plication of topical preparations. Use with wounds | | | | | | of exudate to secondary dressing. May be composed of silicone or the | ulles. |
| SWC-0003 | Non Adherent Sterile Impregnated Gauze 3"x3" | Adaptic | Systagenix 2012 | 3" x 3" 7.6 x 7.6 cm | Y | 4 | | |
| SWC-0076 | Non Adherent Sterile Impregnated Gauze 3"x8" | Adaptic | Systagenix 2015 | 3" x 8" 7.6 x 20.3 cm | Y | 7 | | |
| SWC-0160 | Non Adherent Sterile Silicone Dressing 3"x4" | Mepitel | Mölnlycke 290710 | 3"x 4" 7.5 x 10 cm | Y | 3 | | |
| SWC-0086 | Non Adherent Sterile Impregnated Gauze 4"x4" | Jelonet | Smith & Nephew 7404 | 4"x4" 10 x 10 cm | Y | 4 | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes |
|---------------------|---|------------|------------------------|--|---------|---------------|--|
| | ING DEVICES: Offloading devices are used in ny authorized prescriber (Nurse Practicioner, NSW) | | | | l Conta | act Casts | (TCC), custom braces and ankle and/or foot orthoses. These products can be |
| SWC-0165 | Diabetic Walker Boot System - Small | AirCast XP | DJO Canada 01PD-S | Shoe Size (S) Male 4-7 Female 5-8 | Y | 1 | - The XP Diabetic Walker Boot System incorporates an optimized rocker sole, semi-rigid shells, an |
| SWC-0164 | Diabetic Walker Boot System - Medium | AirCast XP | DJO Canada 01PD-M | Shoe Size (M): Male 7-10 Female 8-11 | Y | 1 | innovative dual density insole and four aircells for maximum off -loading and "total contact" fit of the diabetic foot. Two extra insoles (including the dual- |
| SWC-0163 | Diabetic Walker Boot System - Large | AirCast XP | DJO Canada 01PD-L | Shoe Size (L): Male 10-13 Female 11-15 | Y | 1 | laminate Impax [™] Grid insole), three stockings and one hand bulb with pressure gauge are included. XL available only by special order and may take 3- |
| SWC-0016 | Diabetic Walker Boot System - X Large | AirCast XP | DJO Canada 01PD-XL | Shoe Size (XL) Male 13+ Female 15+ | Y | 1 | 5 days for stock |
| SWC-0017 | Diabetic Walker Insole Replacement Small | AirCast XP | DJO Canada 01PI-S | Small | Y | 1 | The initial order of the Diabetic Walker Boot System |
| SWC-0027 | Diabetic Walker Insole Replacement Medium | AirCast XP | DJO Canada 01PI-M | Medium | Y | 1 | For use with the AirCast XP Diabetic Walker Boot |
| SWC-0041 | Diabetic Walker Insole Replacement Large | AirCast XP | DJO Canada 01PI-L | Large | Y | 1 | XL size available only by special order and may take |
| SWC-0049 | Diabetic Walker Insole Replacement X-Large | AirCast XP | DJO Canada 01PI-XL | Extra Large | Y | 1 | 3-5 days for stock |
| SWC-0162 | Replacement Inflation Device for Diabetic Walker Boot | AirCast XP | DJO Canada 01PG | Each | Y | 1 | Replacement hand bulb with pressure gauge. A single hand bulb with pressure gauge is already included with each order of a Diabetic Walker Boot System. |
| SWC-0166 | Total Contact Cast Kit | Cutimed | BSN Medical 7349900 | 1 Kit | Y | 1 | Kit includes: Stockinette x1, Adehesive Felt x1, Adhesive Foam x1, Padding Roll x2, Soft Fiberglass Cast Tape Rolls x2, Plaster of Paris Slab x1, Fiberglass Cast Tape Roll x3 |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
|---------------------|---|------------------------|------------------------|-----------------------|-----|----------------------------|--|---|
| OTHER & | SPECIALTY: | | | | | | | |
| SWC-0161 | Hydrophilic Wound Paste | Triad | Coloplast 1964 | 71 g Tube | Y | 1 Tube every 2 weeks | Primary mode of action to maintain a moist wound environment, wound debridement and absorbing wound exudate. Sterile and composed zinc oxide paste | And |
| SDR-0011 | Sterile Single Use Scissors | MED-Rx | CHS 85-4202 | 4 1/2" | Y | 7 | Sharp/Blunt. For single-use, sterile techniques only. Not to be reused . Clinic may order 1 / Patient / Visit In-Home may order 1 / Treatment Visit | 20 |
| SDR-0036 | Non-Sterile Reusable Scissors | Fisher Scientific | ALMP12 | 4 1/2" | Y | 1 per Patient Admit | Non-sterile, Single patient use. For non-invasive, non-critical use. May be re-used by the same patient after cleaning or low-level disinfection as per nursing guidelines. Not for use in nursing clinics - Use SDR-0011 Sterile Single Use Scissors instead | Q |
| SDR-0012 | Sterile Plastic Thumb Forceps | MED-Rx | CHS 85-4200 | Each | Y | 1 / Visit | Sterile Single Use only 1 / Patient / Visit | |
| SDR-0025 | Skin Staple Remover | Medline Canada | Medline DYNJ04058 | Each | Y | 1 | Single use; Sterile; Tweezer Style | |
| SDR-0062 | Wound irrigation tip for syringe | MED-RX | MED-RX 20-1104 | 1 Tip | Y | 7 | Single use. Syringe must be ordered separately | |
| SDR-0080 | Tissue Forceps Metal Sterile Disposable | MED-Rx | CHS 85-9000 | 1 x 2 teeth 4-3/4" | Y | 1/Visit | For sharp/surgical debridement and should only be performed by clinicians with this expertise. | |
| SDR-0081 | Curved Iris Scissors Sterile | MED-Rx | CHS 85-9411 | 4 1/2" | Y | 1/Visit | For sharp/surgical debridement and should only be performed by clinicians with this expertise | 90 |
| SWC-0059 | Wound Culture Swab | BD BBL Culture Swab | BD 220105 | Each | Y | 1 | | GREEN: (constant (C-D)) |
| SWC-0064 | Single Use Safety Scalpel - Size 10 | Miltex | Integra 4-510 | Size 10 | Υ | 3 | Sharp/surgical debridement should only be performed by clinicians with this expertise. | |
| SWC-0066 | Flexible Silver Nitrate Caustic 6" Applicator Stick (3) | AMG | AMG Medical 118-395 | Each | Y | 3 | Silver Nitrate 75%, Potassium Nitrate 25% | |
| SWC-0141 | Flexible Silver Nitrate Caustic 6" Applicator Stick (100 pack for Clinic) | AMG | AMG Medical 118-395 | Pack of 100 | N | 1 | For approved ETs nurses and clinic bulk orders only. | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | ΑΡΟ | Weekly Max | Notes | |
|---------------------|--|-------------------|----------------------------|------------------------|----------|---------------|--|--|
| SWC-0121 | No-Sting Barrier Wipe Single Use | Skin-Prep | Smith & Nephew 59420600 | Each 1 mL Wipe | Y | 3 | Helps reduce friction during removal of tapes and films. Can also be used to prepare skin attachment sites for drainage tubes, external catheters, surrounding ostomy sites and adhesive dressings. | Anternational An |
| SWC-0122 | Dermal Curette with 4mm tip | Miltex | Integra 33-54 | Each | Y | 1 | For sharp/surgical debridement and should only be performed by clinicians with this expertise. | 5 |
| SDR-0006 | Cotton Tipped Applicator Plastic Sterile 2/pk | Medline | MDS202095 | 1 Pack of 2 6" Long | Y | 7 | | |
| SDR-0010 | Sterile Tongue Depressor | AMG Medical | 118-421 | Adult | Y | 3 | Adult size Individually wrapped, single patient use. | |
| | ARE KITS: Prepared kits of commonly used ite ems may also be ordered as needed. | ems. Most item | ns in each kit can l | be ordered ind | ividuall | y if an er | tire kit is not required if an individual CHRIS orde | r code is listed below. |
| SDR-0024 | Suture Removal Kit | Medline Canada | Medline DYNDCA1002A | Each | Y | 1 | Single patient use; Sterile, includes Wire Forceps, Gauze, Sponge and Scissors | |
| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
| SMK-0042 Ir | nitial Wound Care Kit: | N/A | Various | Each Kit | Y | 1 | Maximum Quantity: 1 Kit/Patient Admission | |
| SDR-0039 | Non Allergenic Self-Adhesive Dressing Retention | Tape 2" | 1 | 1 | | | Start up kit only. Visiting Nurse to identify and order regular wound care supplies after initial visit | |
| SDR-0064 | Sponge Synthetic Sterile Non-Woven 2"x2" (2 per | Pack) | | 4 | | | | |
| SDR-0065 | Sponge Synthetic Sterile Non-Woven 4"x4" (2 per | Pack) | | 4 | | | No need to order for patients receiving first visit in ambulatory nursing clinics. Clinics are expected to | NO IMAGE AVAILABLE |
| SDR-0006 | Cotton Tipped Applicator Plastic Sterile 2/pk | | | 1 | 1 | | have on hand and replenish stock as needed. | |
| SIV-0084 | Syringe Luer Lok 30 mL | | | 1 | 1 | | | |
| SDR-0062 | Wound irrigation tip for syringe | | | 1 | 1 | | | |
| SWC-0107 | Absorbent Dressing 4"x6" | | | 2 | 1 | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
|---------------------|--------------------------------|------------|----------------------------|----------|-----|---------------|--|-------|
| SWC-0123 Minor | Laceration Tray Kit - Sterile: | Busse | Cardinal Health 312-751 | Each Kit | Y | 1 | Maximum Quantity: 1 Kit/Week | |
| Plastic Tra | у | | | 1 | | | Disposable single use For sharp/surgical debridement and should only be | 0 |
| Gauze Spo | nges | | | 5 | | | performed by clinicians with this expertise. | qri |
| Waterproof | f Sterile Field | | | 1 | | | | |
| Fenestrate | d Waterproof Drape | | | 1 | | | | |
| Absorbent | Towels | | | 2 | | | | |
| Tissue For | ceps | | | 1 |] | | | 60 MR |
| Needle Hol | der | | | 1 | | | | |
| Scissors | | | | 1 | | | | |

NEGATIVE PRESSURE WOUND THERAPY (NPWT)

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes |
|---------------------|-------------|------------|----------------------|------|-----|---------------|-------|
|---------------------|-------------|------------|----------------------|------|-----|---------------|-------|

Negative Pressure Wound Therapy (NPWT): A therapeutic technique indicated for the application of continual or intermittent negative pressure wound therapy to an acute or chronic wound. The device may promote wound healing by the removal of fluids, including wound exudates, irrigation fluids, body fluids and infections materials. The therapy involves the controlled application of sub-atmospheric pressure to the local wound environment, using a sealed wound dressing connected to a vacuum pump. Equipment and supplies are time-limited as per approved policy.

DISPOSABLE SINGLE-USE NPWT: Cardinal NPWT PRO is the primary disposable single-use NPWT device-

Cardinal NPWT PRO Unit, Dressings and Accessories

| SNP-0070 | NPWT PRO Unit | NWPT PRO | Cardinal Health 6708888 | Each | Y | 1 | One unit provided every 42 days for a maximum of 2 per patient admission. Power supply included | |
|----------|--|-----------------------|----------------------------|----------------------------|---|---|--|-----------------------|
| SNP-0072 | NPWT PRO 300mL Canister with Gel | NWPT PRO | Cardinal Health 47-4000 | 300 mL | Υ | 1 | The default primary canister to be used | |
| SNP-0073 | NPWT PRO 500mL Canister with Gel | NWPT PRO | Cardinal Health 47-4500 | 500 mL | Y | 1 | Only order if 300 mL is not enough | |
| SNP-0075 | NPWT PRO Black Foam Dressing Kit - Small | NWPT PRO | Cardinal Health 47-1702 | 10x8x3 cm | Y | 3 | q72 hours. Includes 1 Small Foam Dressing; 1 Drape; 1 Transfer pad with tubing | |
| SNP-0076 | NPWT PRO Black Foam Dressing Kit - Medium | NWPT PRO | Cardinal Health 47-1701 | 20 x 12.5 x 3 cm | Y | 3 | q72 hours. Includes 1 Medium Foam Dressing; 2 Drapes; 1 Transfer pad with tubing | |
| SNP-0077 | NPWT PRO Black Foam Dressing Kit - Large | NWPT PRO | Cardinal Health 47-1700 | 25 x 15 x 3 cm | Y | 3 | q72 hours. Includes 1 Large Foam Dressing; 2 Drapes; 1 Transfer pad with tubing | and the second second |
| SNP-0078 | NPWT PRO Black Foam Dressing Kit - XL | NWPT PRO | Cardinal Health 47-1703 | 58.5 x 33 x 3 cm | Υ | 3 | q72 hours. Includes 1 XL Foam Dressing; 2 Drapes; 1 Transfer pad with tubing | |
| SNP-0079 | NPWT PRO Y-Connector | NWPT PRO | Cardinal Health 47-2500 | Each | Υ | 1 | | 1111111111 |
| SNP-0081 | NPWT PRO White Foam Dressing | NWPT PRO | Cardinal Health 47-1751 | 9.5 × 9.5 x .05 cm | Y | 3 | For sensitive wounds/burns/grafts/tunnelling | |
| SNP-0082 | NPWT PRO Polyurethane Drape | NPWT PRO | Cardinal Health 47-7000 | 25.5 x 33 cm | Y | 2 | Regular drape in stock - Same as provided in NPWT Dressing Kits | |
| SNP-0074 | NPWT PRO SensiSkin Drape | NWPT PRO SensiSkin | Cardinal Health 47-7100 | 14" x 10" 35 cm x 25 cm | Y | 2 | For patients with very sensitive, fragile skin. Hypoallergenic and can help prevent irritation or maceration. Special Order Only - May take up to 5 days for stock. | |
| SNP-0071 | NPWT PRO AC Power Supply | NWPT PRO | Cardinal Health 47-9100 | Each | Y | | Only provided if patient loses original power supply provided with unit (1 per patient admission if required) | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | ΑΡΟ | Weekly Max | Notes | | |
|--------------------------------------|--|------------|----------------------------|--------------------|-----|---------------|---|------------------------|--|
| PICO Negative Pressure Dressing Kits | | | | | | | | | |
| SNP-0069 | PICO-7 Negative Pressure Dressing Kit 4"x8" | PICO 7 | Smith & Nephew 66802002 | 4"x8" 10x20cm | Y | 1 | 2 Dressings in each kit | | |
| SNP-0083 | PICO-7 Negative Pressure Dressing Kit 4"x12" | PICO 7 | Smith & Nephew 66802003 | 4"x12" 10x30 cm | Y | 1 | 2 Dressings in each kit | | |
| SNP-0084 | PICO-7 Negative Pressure Dressing Kit 6"x8" | PICO 7 | Smith & Nephew 66802006 | 6"x8" 15x20 cm | Y | 1 | 2 Dressings in each kit | | |
| SNP-0085 | PICO-7 Negative Pressure Dressing kit 6"x6" | PICO 7 | Smith & Nephew 66802005 | 6"x6" 15x15 cm | Y | 1 | 2 Dressings in each kit | | |
| | E FOAM NPWT: The ActiVAC NPWT sopriate as authorized by the designated Nu | - | | | | | unity care. ActiVAC is only to be used if the options for disp nment with approved policy. | osable single-use NPWT | |
| ETM-1014 | Reusable Foam NPWT LHIN-Owned Unit | ActiV.A.C. | KCI 340013 | Each | Y | | Continuous or Intermittent Modes, Intensity: Low, Medium, High Pressure Settings: -25 to -200 mmHg (-3.3 to -26.6 kPa) Rental only allowed if owned unit is not | | |
| ETM-1010 | Reusable Foam NPWT Rental Unit | ActiV.A.C. | KCI 340013 | Each | Y | | available. Rental unit will automatically be subsitiuted with Owned unit when available as appropriate. | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
|---------------------|---|------------------|----------------------|---------------------|--------|---------------|--|--------------------------|
| | Dressing Kits: Includes the dressings, a | adhesive drapes | | Pad. Excess fluid | and po | - | nfectious materials are removed from the wound through the | Sensa TRAC pad and |
| | VAC GranuFoam Standard Dressing Kit Small | V.A.C. | KCI M8275051 | 10x7.5x3.3 cm | Y | 3 | Includes 1 small black foam dressing, 1 adhesive drape,1 SensaTRAC Pad | |
| SNP-0009 | VAC GranuFoam Standard Dressing Kit Medium | V.A.C. | KCI M8275052 | 18x12.5x 3.3 cm | Y | 3 | Includes 1 medium black foam dressing, 2 adhesive drapes,1 SensaTRAC Pad | |
| SNP-0010 | VAC GranuFoam Standard Dressing Kit Large | V.A.C. | KCI M8275053 | 26x15x3.3 cm | Y | 3 | Includes 1 large black foam dressing, 2 adhesive drapes,1 SensaTRAC Pad | \bigcirc |
| | VAC GranuFoam Standard Dressing Kit XL | V.A.C. | KCI M8275065 | 60x30x1.5 cm | Y | 3 | Includes 1 large black foam dressing, 5 adhesive drapes,1 SensaTRAC Pad. SPECIAL ORDER ITEM - May take up to 5 days for stock. | |
| SNP-0017 | VAC GranuFoam Silver Dressing Kit Small | V.A.C. | KCI M8275098 | 10x7.5x3.3 cm | Y | 3 | Includes 1 small silver foam dressing, 1 adhesive drape,1 SensaTRAC Pad | 1 MI MI MI MI MI MI MI M |
| SNP-0018 | VAC Granufoam Silver Dressing Kit Medium | V.A.C. | KCI M8275096 | 18x12.5x3.3 cm | Y | 3 | Includes 1 medium silver foam dressing, 2 adhesive drapes,1 SensaTRAC Pad | 20/20/2 |
| SNP-0019 | VAC Granufoam Silver Dressing Kit Large | V.A.C. | KCI M8275099 | 26x15x3.3 cm | Y | 3 | Includes 1 large silver foam dressing, 2 adhesive drapes,1 SensaTRAC Pad | Solar Carl |
| SNP-0020 | VAC WhiteFoam Dressing Kit Small | V.A.C. | KCI M8275068 | 10x7.5x1 cm | Y | 3 | Less adherent material helps promote graft take. Includes 1 small PVA foam dressing, 1 adhesive drape, 1 SensaT.R.A.C. Pad | |
| SNP-0022 | VAC WhiteFoam Dressing Kit Large | V.A.C. | KCI M8275067 | 10x15x1 cm | Y | 3 | Less adherent material helps promote graft take. Includes 1 large PVA foam dressing,1 adhesive drape, 1 SensaT.R.A.C. Pad | |
| SNP-0046 | VAC Simplace Dressing Kit Small | V.A.C. | KCI M8275046 | Each | Y | 3 | Includes 1 small Spiral-cut VAC GranuFoam Dressing, 2 adhesive drapes, 1 SensaT.R.A.C. Pad | |
| SNP-0048 | VAC Simplace Dressing Kit Medium | V.A.C. | KCI M8275045 | Each | Y | 3 | Includes 1 medium Spiral-cut VAC GranuFoam Dressing, 2 adhesive drapes, 1 SensaT.R.A.C. Pad | |
| ActiVAC A | Accessories: Additional dsposable comp | oonents to optim | ise individual trea | tment to or when tr | eating | multiple w | vounds. | |
| SNP-0001 | VAC Adhesive Drape | V.A.C. | KCI M6275009 | 30.5 x 26 cm | Y | 1 | Included in the Standard Dressing Kits below. This is an adhesive layer spread across the foam dressing, ensuring a closed, moist wound healing environment. | KCI HKCI |

Each

Υ

3

Y-connector to connect two or more dressings

KCI

M6275066

V.A.C.

SNP-0007

T.R.A.C. Y Connector

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Max | Notes | |
|---------------------|--------------------------------------|------------|----------------------|-----------------|-----|---------------|---|-------------|
| SNP-0058 | VAC WhiteFoam Small (Foam Only) | V.A.C. | KCI M6275033 | 10 x 7.5 x 1 cm | Y | 3 | 1 polyvinyl alcohol (PVA) foam dressing (10 x 7.5 x 1 cm) | |
| SNP-0059 | VAC WhiteFoam Large (Foam Only) | V.A.C. | KCI M6275034 | 10 x 15 x 1 cm | Y | 3 | 1 polyvinyl alcohol (PVA) foam dressing (10 x 15 x 1 cm) | |
| SNP-0052 | 300 mL ActiVAC Canister with Gel | V.A.C. | KCI M8275058 | Each | Y | 1 | 1 canister with gel, tubing, clamp and connector; May be authorized for more than 1/week if the canister is full. | |
| SNP-0060 | VAC SensaT.R.A.C. Pad | V.A.C. | KCI M8275057 | Each | Y | 1 | Includes 1 SensaT.R.A.C. Pad with tubing, clamp and connector | |
| SNP-0061 | T.R.A.C Tubing Cap | V.A.C. | KCI M6275069 | Each | Y | 1 | Cap to secure end of tubing when not in use | E |
| SNP-0067 | VAC GranuFoam Bridge Dressing Kit | V.A.C. | KCI M8275042 | Each | Y | 3 | Each kit includes black foam, bridge dressing; sensa | |
| SNP-0068 | VAC GranuFoam Bridge XG Dressing Kit | V.A.C. | KCI M8275044 | Each | Y | 3 | TRAC pad and adhesive drape. SPECIAL ORDER ITEM - May take up to 5 days for stock. | C Alexandre |

CLEANSERS & ANTISEPTIC SOLUTIONS

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | ΑΡΟ | Weekly Max | Notes |
|---------------------|--|------------|-----------------------------|-----------------|-----|------------|---|
| Cleansers | & Antiseptic Solutions | | | | | | |
| SSO-0001 | Chlorhexidine 2% Alcohol 70% 500 mL Bottle | SoluPrep | 3M 103.15 | 500 mL | Y | 1 | Reclosable bottle. Document date opened on bottle and discard unused product after 30 days |
| SSO-0002 | Hydrogen Peroxide 3% 500mL Bottle | Swan | 189-1000035784 | 500 mL | Υ | 1 | Reclosable bottle. Document date opened on bottle and discard unused product after 7 days |
| SSO-0019 | Normal Saline 500 mL bottle | Baxter | Baxter JF7123 | 500 mL | Y | 1 | When used for for wound care these items are approved only for: |
| SSO-0023 | Normal Saline 118 mL Irrigation Bottle Sterile | Derma | Derma Sciences DUP48501N | 118 mL | Y | 7 | 1) Immunocomprimised patient (i.e. BMT, Oncology, Organ Transplant patient); |
| SSO-0031 | Sterile Water 500 mL bottle | Baxter | Baxter JF7113 | 500 mL | Y | 1 | Well water (untreated) or Environmental considerations (i.e. concerns of hygiene) |
| SSO-0069 | Chlorhexidine 2% Alcohol 70% Wipe | SoluPrep | 3M 101.06 | Each | Y | 28 | Single use antiseptic skin preparation wipe for swabbing sites and ports for vascular access catheters. Single use |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | SoluPrep | 3M 102.03 | Small 1.6 mL | Υ | 8 | Antiseptic skin preparation swabsticks for maintenance of indwelling vascular access catheters. Single use |
| SSO-0052 | Alcohol 70% Wipe | Cardinal | Loris | 30 x 65 mm | Y | 28 | Single use antiseptic topical skin cleanser prior to an injection Medium-sized wipe 70% Isopropyl Alcohol; Single use |
| SSO-0072 | Alcohol 70% Wipe - Box of 200 | Curumu | LP103-03 | | Ν | | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie Sudbury - 3 boxes maximum per month |
| SSO-0054 | Povidone Iodine 10% Wipe | Cardinal | Loris LP109-01 | 30 x 65 mm | Y | | Single use 10% PVP-I antiseptic/germicide prep used for wound clensing or minor invasive procedures, scrapes, cuts or burns. Non- |
| SSO-0067 | Povidone Iodine 10% Swabstick (Single) | Cardinal | Loris LP109-05 | 13 x 25 mm | Y | | irritating and virtually non-stinging for better patient comfort. |
| SSO-0004 | Povidone Iodine 10% 115mL Bottle | Loris | Loris LP109-09 | 115 mL | Υ | 2 | Review manufacturer's product information regarding contraindications. |
| SSO-0011 | Hygeol 1:20 500 mL | N/A | N/A | 500 mL | Y | | Maximum 14/patient |
| SSO-0073 | Wound Irrigation Gel Antimicrobial 30 mL Bottle | Prontosan | B. Braun BB400510 | 30 mL | Y | | For cleansing, moistening and decolonization of acute and chronic skin wounds, 1st and 2nd degree burns. (NERDS: Non-healing, |
| 550-0074 | Wound Irrigation Solution Antimicrobial 350 mL Bottle | Prontosan | B, Braun BB400431 | 350 mL | Υ | | Exudate, Red/bleeding, Debris, Smell). This product is not to be used to replace potable water cleansing |
| Miscellane | eous | | | | | | |
| SOT-0083 | Oral Swab sticks Latex-free Plain Untreated (Box 250) | Toothette | Sage Products 5602UT | Box of 250 | Ν | | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| SOT-0084 | Oral Swab sticks Latex-free Plain Untreated (10 Pack) | Toothette | Sage Products 6072 | Pack of 10 | Y | | Palliative patients only at home/in community |

I.V. SUPPLIES

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|--|--------------|------------------------------|------------------|-----|----------------------|---|
| NEEDLES 8 | & SYRINGES | | | | | | |
| SIV-0079 | Syringe Luer Lok 1mL | | BD | | Y | 7 | Luer Lok tip - Graduation 1/100 mL |
| SIV-0238 | Syringe Luer Lok 1mL - Box of 100 | Luer-Lok | 309628 | 1 mL | N | | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| SIV-0080 | Syringe Luer Lok 3 mL | | BD | | Y | 7 | Luer Lok tip - 0.2 mL graduations |
| SIV-0239 | Syringe Luer Lok 3 mL - Box of 200 | Luer-Lok | 309657 | 3 mL | Ν | 2 Boxes per Month | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| SIV-0081 | Syringe Luer Lok 5 mL | Luer-Lok | BD 309646 | 5 mL | Y | 7 | Luer Lok tip - Graduation 1/5 mL |
| SIV-0082 | Syringe Luer Lock 10 mL | Luer-Lok | BD 302995 | 10 mL | Y | 7 | Luer Lok tip - Graduation 1/5 mL |
| SIV-0084 | Syringe Luer Lok 30 mL | Luer-Lok | BD 302832 | 30 mL | Y | See Notes > | Luer Lok tip - Graduation: 1 mL For IV: 7 / Patient / Week For Wound Care : 1 / Patient / Week In-Home or 1 / Visit in Clinic Also order Wound Irrigation Tip for syringe for each dressing change (SDR-0062) |
| SIV-0085 | Syringe Luer Lok 50 mL | Luer-Lok | BD 309653 | 50 mL | Y | 7 | Luer Lok tip - Graduation: 1 mL , 2 oz in 1/4 oz |
| SIV-0117 | Blunt Fill Needle 18g 1.5" | | BD | 18g 1.5" | Υ | 3 | |
| SIV-0240 | Blunt Fill Needle 18g 1.5" - Box of 100 | BD | 305180 | 1.2 x 40 mm | N | 2 Boxes per Month | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| SIV-0173 | Needle with Syringe 3 mL 23g 1" | Eclipse | BD 305782 | 3mL 23g 1" | Y | 7 | Eclipse Needle with Luer-Lok |
| SIV-0175 | Needle with Syringe 1mL 25g 5/8" | Eclipse | BD 305780 | 1mL 25g 5/8" | Y | 7 | Eclipse Needle with Luer-Lok |
| SIV-0177 | Needle with Syringe 1mL 27g 0.5" | | BD | | Υ | 7 | Eclipse Needle with Luer-Lok |
| SIV-0241 | Needle with Syringe 1mL 27g 0.5" - Box of 50 | Eclipse | 305789 | 1mL 27g 0.5" | N | 3 Boxes per Month | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| SIV-0181 | Needle with Syringe 3mL 25g 5/8" | Eclipse | BD 305781 | 3mL 25g 5/8" | Y | 7 | Eclipse Needle with Luer-Lok |
| SIV-0192 | Gripper Plus 22g x 3/4" Non Y-site | Gripper Plus | Smiths Medical 21-2766-24 | 22 gauge x 3/4" | Υ | 1 | Non-coring needle without Y-site Used to access implanted vascular access devices. |
| SIV-0193 | Gripper Plus 19g x 1.25" Non Y-site | Gripper Plus | Smiths Medical 21-2765-24 | 19 gauge x 1.25" | Υ | 1 | Non-coring needle without Y-site Used to access implanted vascular access devices. |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|--|--------------------------|----------------------|----------------------|-----|-------------------|---|
| SIV-0211 | Blunt Fill Needle with Filter 18g x 1.5" | BD | BD | 18g 1.5" | Υ | 1/Kit | 5 micron filter with blunt fill tip. Only to be used in Symptom Relief Kits for use when withdrawing fluid/medication from an ampoule |
| | Blunt Fill Needle with Filter 18g x 1.5" - Box of 100 | 55 | 305211 | 1.2 x 40 mm | Ν | Month | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| | rth East shall not fund or provide any other co | ncentrations of l | | | | | er's recommendations. Home and Community Care Support |
| SIV-0196 | Heparin Lock Syringe 10 USP units/mL - 5 mL Fill in 10 mL Diameter Syringe | | BD 306414 | 10 USP/mL 5 mL/mL | Υ | 3 | |
| SIV-0199 | Heparin Lock Syringe 10 USP units/mL - 3 mL Fill in 10 mL Diameter Syringe | PosiFlush | BD 306413 | 10 USP/mL 3 mL | Υ | 3 | |
| SIV-0197 | Fill in 10 mL Diameter Syringe | PosiFlush | BD 306424 | 100 USP/mL 5 mL | Υ | 3 | |
| SIV-0200 | Heparin Lock Syringe 100 USP units/mL - 3 mL Fill in 10 mL Diameter Syringe | PosiFlush | BD 306423 | 100 USP/mL 3 mL | Υ | 3 | |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | PosiFlush | BD 306575 | 10 mL | Υ | 14 | |
| I.V. SOLUTI | ONS | | | | | | |
| SSO-0024 | Sodium Chloride 0.9% INJ USP 100 mL Bag | Viaflex | Baxter JB1309M | 100 mL | Υ | 7 | |
| SSO-0025 | Sodium Chloride 0.9% INJ USP 250 mL Bag | Viaflex | Baxter JB1322P | 250 mL | Υ | 7 | |
| SSO-0026 | Sodium Chloride 0.9% INJ USP 500 mL Bag | Viaflex | Baxter JB1323 | 500 mL | Υ | 7 | |
| SSO-0027 | , | Viaflex | Baxter JB1324 | 1000 mL | Υ | 7 | |
| SSO-0029 | | Hospira or equivalent | 954354 | 10 mL | Υ | 7 | Indicated only for diluting or dissolving medications for use in a nebulizer |
| SSO-0032 | Sterlie Water 10 mL Viai | Hospira or equivalent | 954469 | 10 mL | Υ | 7 | Indicated only for diluting or dissolving medications for use in a nebulizer |
| SSO-0070 | Dextrose 5% & Sodium Chloride 0.45% INJ USP 1000 mL Bag | Viaflex | Baxter JB1074 | 1000 mL | Υ | 7 | |
| SSO-0034 | Dextrose 3.3% & Sodium Chloride 0.3% INJ USP 1000 mL Bag | Viaflex | Baxter JB1034 | 1000 mL | Υ | 7 | |
| SSO-0041 | Dextrose 5% INJ USP 1000 mL Bag | Viaflex | Baxter JB0064 | 1000 mL | Υ | 7 | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|---------------------------------------|--------------------------|------------------------------|---|-----|-------------------|--|
| SSO-0071 | 10% Dextrose INJ USP 1000 mL Bag | Viaflex | Baxter JB0164 | 1000 mL | Y | 7 | |
| SSO-0046 | Lactated Ringer's INJ USP 1000 mL Bag | Viaflex | Baxter JB2324 | 1000 mL | Y | 7 | |
| I.V. CATHET | ERS & ACCESS | | | - | - | - | |
| SIV-0150 | IV Catheter 22g Blue | Nexiva | BD 383512 | 22 gauge1.00" 0.9 mm x 25 mm | Y | 3 | Closed IV catheter system-Single port |
| SIV-0151 | IV Catheter 24g Yellow | Nexiva | BD 383511 | 24 gauge 0.75" 0.7 mm x 19 mm | Y | 3 | Closed IV catheter system-Single port |
| SIV-0228 | IV Catheter 18g Green | Nexiva | BD 383519 | 18 gauge 1.25" 1.1 mm x 32 mm | Ν | 2 | Closed IV catheter system-Single port For MAiD Patients as approved by Care Coordinator |
| SIV-0229 | IV Catheter 20g Pink | Nexiva | BD 383517 | 20 gauge 1.25" 1.3 mm x 32 mm | Ν | 4 | Closed IV catheter system-Single port For MAiD Patients as approved by Care Coordinator |
| SIV-0187 | IV Catheter Set - 22g | Saf-T-Intima | BD 383329 | 22 gauge 0.75" 0.9 mm x 19 mm | Y | 3 | Safety system with wings, extension tubing and Y adapter |
| SIV-0203 | IV Catheter Set - 24g | 0 (| BD | 24 gauge 0.75" | Y | 3 | Safety system with wings, extension tubing and Y adapter |
| SIV-0243 | IV Catheter Set - 24g - Box of 25 | Saf-T-Intima | 383319 | 0.9 mm x 19 mm | Ν | | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| I.V. GRAVIT | Y SETS AND EXTENSION SETS | | | | | | |
| SIV-0225 | Continuous Flow IV Solution Set | Clearlink Continu-Flo | Baxter 2C8519 | 10 drops / mL 112" (2.8 m) | Y | 7 | 2 Luer Activated Valves, Male Luer Lock Adapter, Approximate Volume 10 drops per mL |
| SIV-0226 | Secondary Medication Set | Clearlink | Baxter JC7453 | 10 drops / mL 37" (94 cm) | Y | 7 | With Male Luer lock adapter and hanger |
| SIV-0009 | Needle-Free Connector Valve | MaxZero | Cardinal Health MZ1000-07 | Length: 3 cm Priming Volume 0.19 mL | Y | 7 | Clear needleless connector |
| SIV-0012 | Catheter Extension Set | MaxZero | Cardinal Health MZ5301 | Length: 18 cm Priming Volume 0.3 mL | Y | | Extension set with clear needle-free connector valve |
| SIV-0224 | Y Connector IV Extension Set | MaxZero | Cardinal Health MZ5307 | 7" / 8 cm PV: 0.8 mL | Y | 1 | Bifuse pressure rated extension set, 2 bonded MaxZero connectors, 2 slide clamps, spin male luer lock |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|---|-------------------------|------------------------|----------------------------------|--------|----------------------|---|
| | REMENT: Securement devices and dressings u rm IV Advanced and IV3000 is only to be autho | | | | | | ntral catheter insertion sites and wounds. The primary product is y not appropriate due to sensitivity. |
| SDR-0094 | Tegaderm IV Advanced Securement Dressing | | 3M | 2.5 "x 2.75" | Y | 3 | For peripheral IVs |
| SIV-0244 | Tegaderm IV Advanced Securement Dressing - Box of 100 | Tegaderm IV | 1683 | 6.5 cm x 7 cm | N | 2 Boxes per Month | For Residential Hospice bulk orders only in North Bay, Sault Ste. Marie, Sudbury - 3 boxes maximum per month |
| SIV-0231 | PICC/CVC Securement Device and Dressing Kit (Regular) | Tegaderm IV Advanced | 3M 1837-2100 | 3 .5" x 4.5" 8.5 cm x 11.5 cm | Y | 1 | To secure Peripherally Inserted Central Catheters (PICC's) and short term Central Venous Catheters (CVC's) Designed to be worn for up to 7 days |
| SIV-0232 | PICC/CVC Securement Device and Dressing Kit (Large) | Tegaderm IV Advanced | 3M 1839-2100 | 4" x 6" 10 cm x 15.5 cm | | | Removes cleanly, without causing patients undue pain or distress Accommodates the majority of PICC and CVC catheters up to & incl. 12 FR |
| SIV-0236 | IV3000 Securement Dressing | IV 3000 1Hand | Smith & Nephew 4008 | 4" x 4-3/4" 10cm x 12cm | Y | 2 | Only to be ordered if Tegaderm IV dressings are not appropriate. Dressing only for Central/PICC/midline catheters. |
| SIV-0134 | StatLock IV Catheter Securement Device | StatLock PICC Plus | Bard PIC0220 | Each | Y | 1 | Secures catheter in place when using the IV3000 dressing <u>when required</u> to secure PICC line. Accommodates a wide variety of PICC Catheters. Tricot Anchor Pad, Sliding Post, Sterile. |
| NOTE: Addi | itional items are also available in the TRACHEOT | OMY / DRAIN S | PONGES in the Dre | essings, Bandages | , Tape | , Compress | \mathbf{H} |
| INFUSION | PUMPS & SUPPLIES | | | | | | |
| ETM-1005 | IV pole - portable | | Rental | equipment | | | |
| ETM-1009 | Ambulatory infusion pump (includes pump, power supply, small carrying pouch and batteries) | | Rental | equipment | | | All service areas |
| ETM-1012 | Large carry case - for ambulatory infusion pump (if required)(For infusions requiring 250 mL or 500 mL bags) | | Rental | equipment | | | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| ETM-1016 | Remote Dose Extension Cord for ETM-1009 IV Pump | | Rental | equipment | | | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| ETM-0028 | Sapphire IV Pump Mini-Cradle Mount for Pole | | Rental | equipment | | | Cochrane sub-region/Timmins Branch area (Bayshore Specialty Rx) |
| ETM-0030 | Sapphire IV Pump Wearable Pouch 500 ml | | Rental | equipment | | | Cochrane sub-region/Timmins Branch area (Bayshore Specialty Rx) |
| ETM-0034 | Sapphire IV Pump Bolus Handle | | Rental | equipment | | | Cochrane sub-region/Timmins Branch area (Bayshore Specialty Rx) |
| SOT-0026 | Replacement batteries for CADD IV pump | | | 4- AA's | Y | 1 | Single use batteries for CADD infusion pump- Solis Not used in Timmins Branch/Cochrane sub-region |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|---|---------------------------|------------------------------|-------------------------|-----|-------------------|--|
| | | Bodyguard | CME BBA120-161XPS | 80" | Υ | 7 | Timmins Branch/Cochrane sub-region (Bayshore Specialty Rx) |
| SIV-0042 | Ambulatory Infusion Pump Administration Set | Sapphire | AP-407 12000-000-0008 | 112" 286 cm | Υ | 7 | |
| | | CADD | Smiths Medical 21-7022-24 | Each | Υ | 7 | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| SIV-0043 | Ambulatory Infusion Pump High Volume Administration Set | CADD | Smiths Medical 21-7345-24 | Each | Y | 3 | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| | | Bodyguard | CME A120-112XSXFK | 110" | Υ | 7 | Timming Branch/Cookrang out region (Doughors Specialty Du) |
| SIV-0044 | Ambulatory Infusion Pump Administration Set Filtered | Sapphire | AP 206 12000-000-0005 | 116" 295 cm | Y | 7 | Timmins Branch/Cochrane sub-region (Bayshore Specialty Rx) |
| | | CADD | Smiths Medical 21-7346-24 | Each | Y | 7 | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| 0045 | Anabulatan Jufusian Duma Estanaian Tuking Cat | ICU | ICU Medical B1327 | 20 inches | - Y | 3 | Timmins Branch/Cochrane sub-region (Bayshore Specialty Rx) |
| SIV-0045 | Ambulatory Infusion Pump Extension Tubing Set | CADD | Smiths Medical 21-7045-24 | 30 Inches | Y | 3 | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| SIV-0046 | Ambulatory Infusion Pump Extension Tubing 45 inches | CADD | Smiths Medical 21-7046-24 | 45 inches | Y | 3 | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| SIV-0047 | Ambulatory Infusion Pump Medication Cassette 50ml | | | 50 mL | Y | 7 | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| SIV-0048 | Ambulatory Infusion Pump Medication Cassette 100ml | | | 100 mL | Υ | 7 | Sudbury, North Bay, Kirkland Lake, Sault Ste. Marie, Parry Sound Branch Areas (Robinsons Pharmacy and Medicine Shoppe) |
| SIV-0202 | IV Medication Refill Solution Bag | | | Each | Υ | 7 | Solution bag - not a cassette. For IV Med Refill orders only |
| physician, el | | vided if medica | tions are funded ur | | | • | ump. If a single-use pump is requested by the patient or 3). Other costs related to elastomeric pumps that are not covered |
| SIV-0233 | Elastomeric Single Use Disposable IV Pump | EasyPump or equivalent | B.Braun | Each | Ν | 1 | SPECIAL ORDER ITEM: Care Coordinator to confirm ODB coverage and timing for availability with the approved IV pharmacy prior ordering pump. The pump type/size is selected by pharmacy per medication type. |
| SIV-0234 | Elastomeric IV Pump Catheter Extension Set | Max Plus Clear | MP9014C | 14"/36 cm PV 0.75 ml | Y | 3 | For use with SIV-0233 Elastomeric pump only. |
| SIV-0235 | Elastomeric IV Pump Belt Carrying Bag | EasyPump or equivalent | B.Braun | Each | Ν | 1 / patient | For use with SIV-0233 Elastomeric pumps only. |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|---|----------------------|------------------------|-------------------------|---------|-------------------|--|
| BLOOD CO | DLLECTION SUPPLIES: For blood collection for | or patients eligit | ole for in-home La | b Services as initi | iated b | y Care Coo | rdinator only as per policy |
| SIV-0137 | Stretch Tourniquet | Vacutainer | BD 367203 | 18" x 1" | Y | 1 | |
| SIV-0145 | Blood collection infusion set with luer adapter 23g | Safety-Lok | BD 367283 | 23 g x 0.75" Needle | Y | 3 | 12" Tubing Manually activated safety shield to fully cover needle. |
| SIV-0204 | Blood collection infusion set with luer adapter 21g | Safety-Lok | BD 367281 | 21 g x 0.75" Needle | Υ | 3 | 12" Tubing Manually activated safety shield to fully cover needle |
| SIV-0212 | Blood Tube holder Single Use Nonstackable | Vacutainer | BD 364815 | Each | Y | 1 | Single patient use only - Discard the entire device. Needle is not part of the tube holder. One Blood Tube Holder can be used for all the blood work per patient |
| SIV-0205 | Venous Blood Collection Tube - Lavender | Vacutainer | BD 367841 | 13 x 75 mm x 2.0 mL | Y | 3 | |
| SIV-0206 | Venous Blood Collection Tube - Yellow | Vacutainer | BD 367986 | 13 x 100 mm x 5.0 mL | Y | 3 | |
| SIV-0207 | Venous Blood Collection Tube - Blue | Vacutainer | BD 363083 | 13 x 75 mm x 2.7 mL | Υ | 3 | |
| SIV-0208 | Venous Blood Collection Tube - Green | Vacutainer | BD 367962 | 13 x 100 mm x 4.5 mL | Y | 3 | |
| SIV-0209 | Venous Blood Collection Tube - Red | Vacutainer | BD 367815 | 13 x 100 mm x 6.0 mL | Y | 3 | |
| SIV-0210 | Lab Specimen Bag | Fisher Scientific | Fisherbranded 22310004 | 6" x 10" | Υ | 1 | Transparent Blue, Non-Sterile Single Use Disposable With separate outside pocket holds the requisition slip |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes | | | | | |
|---------------------|--|------------------------|--------------------------------------|------------------------------|-----|-------------------|--|--|--|--|--|--|
| OTHER I.V. | THER I.V. SUPPLIES & ACCESSORIES | | | | | | | | | | | |
| SIV-0132 | Sterile non-fenestrated drape 18x26" | Busse or equivalent | Cardinal Health 696 or equivalent | 18 x 26" | Y | 2 | Multi-layered, absorbent, all-purpose sterile field | | | | | |
| SIV-0213 | Vacuum Drainage Bottle 1000mL with 120cm Tubing Set | PreVac MedRx | CHS S55148 | 1000 mL | Y | 3 | Sterile bottle includes attached drain line with male luer connector and roller camp (drain at 12 PSI). Not to be used for body fluid samples for diagnostic testing | | | | | |
| SIV-0214 | Drainage Procedure Kit (Excludes Bottle) | MedRx | CHS 2976C | 1 Kit | Y | 3 | Order SIV-0213 drainage bottle separately. Includes dressings, tubing, alcohol wipes, vented luer cap, tube clamp, and chemo-rated Nitrile gloves. Sterile. | | | | | |
| SIV-0215 | Sterile Vented Luer Cap | LSL Health Care | CHS 797A | Each | Y | 3 | Not required if ordering Kit 2976C (SIV-0214) For pts with luer lock connections | | | | | |
| SIV-0245 | Drainage Bag 2000mL with 120cm Tubing Set | MedRx | CHS M04200 | 2000 mL | Y | 7 | Includes 120 cm detachable drain line with rotating male luer fitting; roller clamp; built in rotating hanger; integrated drip chamber and anti-reflux valve | | | | | |
| SIV-0217 | Safety Drain Line Adapter | MEDRx | CHS SA5200 | Each | Y | 3 | For use ONLY if patient has Pleurex™ catheter | | | | | |
| SIV-0227 | Universal Catheter Cap | Medionics | CHS S55300 | Each | Y | 3 | For use ONLY if patient has Pleurex™ catheter | | | | | |
| SIV-0131 | Sharps Container Phlebotomy | SharpSafety | Cardinal Health Z8900SA | 6.25x4.25x4.25" 1 Litre | Y | 1 | Accommodates a variety of Sharps; Temporary closure for security between uses; Lids lock for safe transport and disposal. | | | | | |
| SIV-0218 | Sharps Container Chemotherapy | ChemoSafety | Cardinal Health 8982CRH | 10"x7.25"x10.5" 7.6 Litre | Y | 1 | For the disposal of chemotherapy or other cytotoxic waste or for MAiD patients as approved by the Care Coordinator. Temporary closure and locking lid for safe transport and disposal. | | | | | |
| SIV-0219 | Spill Kit for Cytotoxic Drugs | Safetec | Safetec 14399 | Each | Y | 1 | Protects and cleans up spills of chemo or other hazardoudous medications. All inclusive kits contains all needed supplies and instructions. | | | | | |

| CHRIS Order Code | Description | Brand Name | APO | Weekly Maximum | Notes | |
|---------------------|---|-------------|--|-------------------|-------|--|
| | epared kits of commonly used items. Most iter lso be ordered as needed. | not require | ad and the individual CHRIS order code is listed below. Additional | | | |
| SMK-0029 PI | ICC Dressing Kit: | 1 | Maximum Quantity: 1 Kit/Week | | | |
| SDR-0005 | Dressing Tray Kit | | 1 | | | All PICC/central lines require pump for infusion therapy. A large Ambulatory Infusion Pump carrying bag may be required, if necessary. |
| SDR-0028 | Tubular Net Elastic Bandage - Elbow & Arm | | 1 m | | | Add: 2 of Ambulatory Infusion Pump Administration Sets (SIV-0042) OR |
| SIV-0231 | PICC/CVC Securement Device and Dressing Kit (I | Regular) | 1 | | | Ambulatory Infusion Pump Administration Sets Filtered for administration of Cephalosporins, Vancomycin, Ceftazidine (SIV-0044) and 2 Alcohol |
| SIV-0009 | Needle-Free Connector Valve | | 1 | | | Wipes (SS0-0052) AND 2 Sodium Chloride 0.9% Syringes (SSO-0021) per day of infusion |
| SOT-0001 | Disposable Ear Loop Surgical Face Mask | | 2 | | | Add: Portacath Gripper item SIV-0192 or SIV-0193 as required. |
| SOT-0021 | Sterile Nitrile Gloves Medium (Pair) | | 1 | | | Add: Prefilled Heparin Lock Syringes as per nursing guidelines (See IV, Needles, Syringes) |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | 2 | | | Dressing Tray SDR-0005 may be ordered individually if this entire PICC |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | | 2 | | | kit is not required. |
| For double lu | umen add the following to your order: | | | | | |
| SIV-0009 | Needle-Free Connector Valve | | 1 | Y | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | 1 | Y | | |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | | 2 | Y | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---|---|------------|----------------------|------|----------------|--|---|
| SMK-0039 Subcutaneous Administration of Pain Management Medications via Infusion Pump: Initial Start Up Kit | | | | | | | Maximum Quantity: 1 Kit/Week |
| SWC-0072 | Transparent Film Dressing 10 x 12cm | | | 2 | | | Add: Extension tubing (SIV-0045 or SIV-0046) as needed if being used with medication cassette. |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | |
| SIV-0203 | IV Catheter Set - 24g | | | 2 | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 1 | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 4 | | | |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | | | 1 | | | |
| SSO-0052 | Alcohol 70% Wipe | | | 4 | | | |
| SMK-0040 P | Peripheral IV Start Kit: | | | | Y | 2 | Maximum Quantity: 2 Kits/Week |
| SDR-0003 | Gauze Sterile Woven 2"x2" (2 per Pack) | | 2 | | Set (SIV-0042) | | |
| SDR-0028 | Tubular Net Elastic Bandage - Elbow & Arm | | | 1 m | | Ambulatory Infusion Pump Administration Set Filtered for of Cephalosporins, Vancomycin, Ceftazidine (SIV-0044) <u>AND</u> 2 Chlorhexidine 2% Alcohol 70% Wipes (SS0-0069) <u>AND</u> | <u>OR</u> Ambulatory Infusion Pump Administration Set Filtered for administration |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | AND |
| SIV-0009 | Needle-Free Connector Valve | | | 1 | | | AND |
| SIV-0137 | Stretch Tourniquet | | 1 | | | 2 Sodium Chloride 0.9% Syringes (SSO-0021) <u>per day</u> of infusion | |
| SIV-0150 | IV Catheter 22g Blue | | | 2 | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 4 | | | |
| SSO-0052 | Alcohol 70% Wipe | | | 4 | 1 | | |
| SDR-0094 | Tegaderm IV Advanced Securement Dressing | | | 2 | 1 | | |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | | | 2 | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Maximum | Notes |
|--|---|------------|----------------------|------|-----|---------|---|
| SMK-0041: Gravity Peripheral IV Administration Kit | | | | | | | Maximum Quantity: 2 Kit/Week |
| SIV-0137 | Stretch Tourniquet | | | 1 | | | Provides 3 days of infusion For increased frequency of medication administration beyond once daily, |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | order additional Alcohol wipes (SSO-0052) and Prefilled Sodium Chloride 0.9% Syringes (SSO-0021) per day of infusion |
| SDR-0003 | Gauze Sterile Woven 2"x2" (2 per Pack) | | | 2 | | | |
| SIV-0225 | Continu-Flo solution set | | | 3 | | | If additional line required add SIV-0226 Secondary IV set to the order |
| SSO-0024 | Sodium Chloride 0.9% INJ USP 100 mL Bag | | | 3 | | | Additional catheters may be ordered as needed |
| SIV-0151 | IV Catheter 24g Yellow | | | 2 | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 1 | | | |
| SDR-0028 | Tubular Net Elastic Bandage - Elbow & Arm | | | 1 m | | | |
| SDR-0094 | Tegaderm IV Advanced Securement Dressing | | | 1 | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 4 | | | |
| SSO-0052 | Alcohol 70% Wipe | | | 6 | | | |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | | | 6 | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes | | |
|---|---|------------|----------------------|------|-----|--------------------|--|--|--|
| SYMPTOM RELIEF SUPPLY KITS: Geographic area-specific medical supplies kits for pain and symptom management for SRC 95 End-of-Life patients. Prepared kits of commonly used items. Items in each kit can be ordered individually if an entire kit is not required. Additional items may also be ordered as needed. | | | | | | | | | |
| SIV-0223: Symptom Relief Supply Kit - North Bay | | | | | | | Medication order must be provided to the geographically designated approved Infusion Supplies pharmacy provider with the supplies order. The approved Infusion Supplies provider will deliver medication and | | |
| SDR-0086 | White Surgical Tape 1" Single Use Roll 1 | | | | | | supply kit to the patient's home or alternate treatement address specified in the order | | |
| SDR-0064 | Sponge Synthetic Sterile Non-Woven 2"x2" (2 per Pack) 2 | | | | | | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 2 | | | | | |
| SIV-0117 | Blunt Fill Needle 18g 1.5" | | | 5 | | | | | |
| SIV-0175 | Needle with Syringe 1mL 25g 5/8" | | | 10 | | | | | |
| SIV-0203 | IV Catheter Set - 24g | | | 2 | | | | | |
| SSO-0052 | Alcohol 70% Wipe | | | 20 | | | | | |
| SWC-0072 | Transparent Film Dressing 10x12cm | | | 2 | | | | | |
| SIV-0220: Sy | ymptom Relief Supply Kit - Sudbury | | | | Y | order. Includes he | Medication order must be provided to the pharmacy with the supplies order. Includes home delivery within Greater Sudbury only: Designated approved Infusion Supplies pharmacy provider will deliver | | |
| SIV-0131 | Sharps Container Phlebotomy | | | 1 | | | medication and supply kit to the patient's treatement address as specified in the order. | | |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | | | |
| SDR-0094 | Tegaderm IV Advanced Securement Dressing | | | 4 | | | | | |
| SIV-0203 | IV Catheter Set - 24g | | | 4 | 1 | | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 3 | | | | | |
| SIV-0175 | Needle with Syringe 1mL 25g 5/8" | | | 25 | 1 | | | | |
| SSO-0052 | Alcohol 70% Wipe | | | 20 | 1 | | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|---|------------|----------------------|------|-----|-------------------|--|
| SIV-0221: Sy | SIV-0221: Symptom Relief Supply Kit - Timmins | | | | | | Medication order must be provided to the Designated approved Infusion Supplies pharmacy provider and the Home and Community Care North East Internal Medical Supplies prepares this supply kit |
| SSO-0052 | Alcohol 70% Wipe | | | 25 | | | Last internal medical Supplies prepares this supply kit |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 4 | | | |
| SIV-0079 | Syringe Luer Lok 1mL | | | 25 | | | |
| SIV-0080 | Syringe Luer Lok 3 mL | | | 15 | | | |
| SIV-0117 | Blunt Fill Needle 18g 1.5" | | | 40 | | | |
| SIV-0203 | IV Catheter Set - 24g | | | 4 | | | |
| SWC-0072 | Transparent Film Dressing 10x12cm | | | 4 | | | |
| SIV-0131 | Sharps Container Phlebotomy | | | 1 | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 4 | | | |
| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|---|------------|----------------------|------|-----|-------------------|--|
| SIV-0230: Sy | ymptom Relief Supply Kit - Kirkland Lake | | | | N | 1 | Medication order must be provided to the Designated approved Infusion Supplies pharmacy provider and the Home and Community Care North East Internal Medical Supplies prepares this supply kit |
| SSO-0052 | Alcohol 70% Wipe | | | 25 | | | Last internal medical Supplies prepares this supply fit |
| SIV-0211 | Blunt Fill Needle with Filter 18g x 1.5" | | | 7 | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 4 | | | |
| SDR-0064 | Sponge Synthetic Sterile Non-Woven 2"x2" (2 per | Pack) | | 4 | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 4 | | | |
| SIV-0079 | Syringe Luer Lok 1mL | | | 40 | | | |
| SIV-0175 | Needle with Syringe 1mL 25g 5/8" | | | 20 | | | |
| SIV-0203 | IV Catheter Set - 24g | | | 4 | | | |
| SWC-0072 | Transparent Film Dressing 10x12cm | | | 4 | | | |
| SIV-0131 | Sharps Container Phlebotomy | | | 1 | | | |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|--|------------|----------------------|------|-----|-------------------|--|
| SIV-0237: Sy | ymptom Relief Supply Kit - Manitoulin | | | | Y | | The medication order may be provided by any local pharamcy, but the supplies are provided by the designated approved Infusion Supplies |
| SSO-0052 | Alcohol 70% Wipe | | | 25 | | | pharmacy provider |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 4 | | | |
| SIV-0079 | Syringe Luer Lok 1mL | | | 25 | | | |
| SIV-0080 | Syringe Luer Lok 3 mL | | | 15 | | | |
| SIV-0117 | Blunt Fill Needle 18g 1.5" | | | 40 | | | |
| SIV-0203 | IV Catheter Set - 24g | | | 4 | | | |
| SWC-0072 | Transparent Film Dressing 10x12cm | | | 4 | | | |
| SIV-0131 | Sharps Container Phlebotomy | | | 1 | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 4 | | | |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | |
| SIV-0175 | Needle with Syringe 1mL 25g 5/8" | | | 20 | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|---|------------|----------------------|------|-----|-------------------|--|
| SIV-0222 Sy | mptom Relief Supply Kit - Sault Ste Marie | | | | Y | 1 | The medication order may be provided by any local pharamcy, but the supplies are provided by the designated approved Infusion Supplies |
| SWC-0072 | Transparent Film Dressing 10x12cm | | | 3 | | | pharmacy provider |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | |
| SIV-0203 | IV Catheter Set - 24g | | | 3 | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 3 | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 4 | | | |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | | | 1 | | | |
| SSO-0052 | Alcohol 70% Wipe | | | 30 | | | |
| SIV-0117 | Blunt Fill Needle 18g x 1.5:" | | | 30 | | | |
| SIV-0079 | Syringe Luer Lok 1mL | | | 30 | | | |
| SIV-0080 | Syringe Luer Lok 3 mL | | | 5 | | | |
| SIV-0131 | Sharps container Phelbotomy | | | 1 | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|---|------------|----------------------|------|---|-------------------|---|
| MAIDKIT: MA | AiD Supplies Kit | | N | 1 | Pre-assembled kit of supplies required for use by the physican with eligible, active patients approved for MAiD as per approved Home and Community Care Support Services North East policy. | | |
| SDR-0003 | Gauze Sterile Woven 2"x2" (2 per Pack) | | | 4 | | | |
| SDR-0004 | Gauze Sterile Woven 4"x4" (2 per Pack) | | | 4 | | | When placing an order from Robinson's or Internal Medical Supplies, add the delivery code #FC-1033 for LHIN-funded Home Delivery for |
| SDR-0028 | Tubular Net Elastic Bandage - Elbow & Arm | | | 1m | | | MAiD to this order. This Finance Charge is <u>in addition to</u> the appropriate Urban or Rural Delivery Finance Charge and <u>only applicable to this kit</u> . |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | |
| SDR-0094 | Tegaderm IV Advanced Securement Dressing | | | 2 | | | |
| SIV-0009 | Needle-Free Connector Valve | | | 2 | | | |
| SIV-0012 | Catheter Extension Set 18 cm | | | 2 | | | |
| SIV-0218 | Sharps Container Chemotherapy | | | 1 | | | |
| SIV-0137 | Stretch Tourniquet | | | 1 | | | |
| SIV-0224 | Y Connector IV Extension Set | | | 2 | | | |
| SIV-0228 | IV Catheter 18g Green | | | 2 | | | |
| SIV-0229 | IV Catheter 20g Pink | | | 4 | | | |
| SOT-0096 | Non Sterile Nitrile Gloves Large (Pair) | | | 4 | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 8 | | | |
| SSO-0021 | Sodium Chloride 0.9% 10 mL Pre-Filled Syringe | | | 6 | | | |
| SSO-0052 | Alcohol 70% Wipe | | | 8 | | | |

OSTOMY SUPPLIES & ACCESSORIES

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | ΑΡΟ | Weekly Maximum | Notes |
|---|--|---|--|---|----------|----------------------------------|--|
| Authorization private pay c Ostomy SI | n for Home and Community Care Support Serv options, insurance coverage and ADP applicatio | ices North East ins. Coloplast i flanges (also ca | t-funded supplies is the primary vend lled "baseplates") t | as per approved Ostom dor. hat attach directly to the | skin. Me | y and only for easure the sto | agement of urostomies, ileostomies and colostomies. short-term treatment to allow patient some time to coordinate ma accurately and select a cut-to-fit barrier opening which fits |
| SOS-0068 | 2 Piece Click Barrier Green Flat Cut-to-Fit | SenSura Mio Click | Coloplast 10502 | Coupling/Flange 40 mm (1 9/16") | Y | | Non-Convexity, Stoma size 10-35mm (3/8"- 1 3/8") |
| SOS-0025 | 2 Piece Click Barrier Red Flat Cut-to-Fit | SenSura Mio Click | Coloplast 10512 | Coupling/Flange 50 mm (2") | Y | | Non-Convexity, Stoma size 10-45mm (3/8"- 1 3/4") |
| SOS-0026 | 2 Piece Click Barrier Yellow Flat Cut-to-Fit | SenSura Mio Click | Coloplast 10531 | Coupling/Flange 70 mm (2 3/4") | Y | | Non-Convexity, Stoma size 10-65mm (3/8"- 2 1/5") |
| OS-0027 | 2 Piece Click Barrier Green Convex Light Cut-to- Fit | SenSura Mio Click | Coloplast 16901 | Coupling/Flange 40 mm (1 9/16") | Y | | Light Convexity, Stoma size 10-20mm (3/8"- 3/4") |
| SOS-0028 | 2 Piece Click Barrier Red Convex Light Cut-to-Fit | SenSura Mio Click | Coloplast 16911 | Coupling/Flange 50 mm (2") | Y | | Light Convexity, Stoma size 15-30mm (5/8"- 1 1/8") |
| SOS-0029 | 2 Piece Click Barrier Blue Convex Light Cut-to- Fit | SenSura Mio Click | Coloplast 16921 | Coupling/Flange 60 mm (2 3/8") | Y | 5 / Week | Light Convexity, Stoma size 15-40mm (5/8"- 1 9/16") |
| OS-0030 | 2 Piece Click Barrier Yellow Convex Light Cut-to- Fit | SenSura Mio Click | Coloplast 16931 | Coupling/Flange 70 mm (2 3/4") | Y | Maximum 10 | Light Convexity, Stoma size 15-53mm (5/8"- 2 1/16") |
| OS-0031 | 2 Piece Click Barrier Green Convex Deep Cut-to- Fit | SenSura Mio Click | Coloplast 16941 | Coupling/Flange 40 mm (1 9/16") | Y | per patient admission | Deep Convexity, Stoma size 10-20mm (3/8"- 3/4") |
| OS-0032 | 2 Piece Click Barrier Red Convex Deep Cut-to-Fit | SenSura Mio Click | Coloplast 16951 | Coupling/Flange 50 mm (2") | Y | of any combination | Deep Convexity, Stoma size 15-30mm (5/8"- 1 1/8") |
| SOS-0033 | 2 Piece Click Barrier BlueConvex Deep Cut-to-Fit | SenSura Mio Click | Coloplast 16961 | Coupling/Flange 60 mm (2 3/8") | Y | of sizes | Deep Convexity, Stoma size 15-40mm (5/8"- 1 9/16") |
| OS-0034 | 2 Piece Click Barrier Yellow Convex Deep Cut-to- Fit | SenSura Mio Click | Coloplast 16971 | Coupling/Flange 70 mm (2 3/4") | Y | | Deep Convexity, Stoma size 15-53mm (5/8"- 2 1/16") |
| OS-0017 | 2 Piece Hollister Barrier Green Flat Cut-to-fit | CeraPlus | Hollister 11202 | Coupling/Flange 44 mm (1 3/4") | Y | | Non-Convexity, Stoma size <u>up to</u> 32mm (1 1/4") Infused with ceramide to protect peristomal skin |
| OS-0018 | 2 Piece Hollister Barrier Red Flat Cut-to-fit | CeraPlus | Hollister 11203 | Coupling/Flange 57 mm (2 1/4") | Y | | Non-Convexity, Stoma size <u>up to</u> 44mm (1 3/4") Infused with ceramide to protect peristomal skin |
| OS-0019 | 2 Piece Hollister Barrier Green Convex Cut-to-fit | CeraPlus | Hollister 11402 | Coupling/Flange 44 mm (1 3/4") | Y | | Convex, Stoma size <u>up to</u> 25mm (1") Infused with ceramide to protect peristomal skin |
| SOS-0020 | 2 Piece Hollister Barrier Red Convex Cut-to-fit | CeraPlus | Hollister 11403 | Coupling/Flange 57 mm (2 1/4") | Υ | | Convex, Stoma size <u>up to</u> 38mm (1 1/2") Infused with ceramide to protect peristomal skin |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes | |
|---------------------|--|---------------------------|----------------------|-----------------------------------|----------|-----------------------------------|--|----|
| | Pouches for 2 Piece System: Order with a drain valve at the bottom to empty the bag. | appropriate bra | nd and size Oston | ny Skin Barrier A draina | able pou | ch that conne | cts to the barrier flange to collect urine with | - |
| SOS-0035 | 2 Piece SenSura Mio Click Urostomy Pouch Green | SenSura Mio Click Maxi | Coloplast 11492 | Coupling/Flange 40mm (1 9/16") | Y | | Bag length 26cm (10 1/4"), | O. |
| SOS-0036 | 2 Piece SenSura Mio Click Urostomy Pouch Red | SenSura Mio Click Maxi | Coloplast 11498 | Coupling/Flange 50mm (2") | Υ | Maximum 10 per patient | Transparent | |
| SOS-0021 | 2 Piece Hollister Urostomy Pouch Green | Hollister | Hollister 18902 | Coupling/Flange 44mm (1 3/4") | Y | admission of any | | R |
| SOS-0022 | 2 Piece Hollister Urostomy Pouch Red | Hollister | Hollister 18903 | Coupling/Flange 57mm (2 1/4") | Υ | combination of sizes | | - |
| | Pouches for 2 Piece System: Order with a to the barrier flange. An integrated closure or clar | | | | | ch with a large | r opening for ileostomies and colostomies | - |
| SOS-0037 | 2 Piece SenSura Mio Click Drainable Pouch Green | SenSura Mio Click Maxi | Coloplast 11463 | Coupling/Flange 40mm (1 9/16") | Y | | | Ö |
| SOS-0038 | 2 Piece SenSura Mio Click Drainable Pouch Red | SenSura Mio Click Maxi | Coloplast 11473 | Coupling/Flange 50mm (2") | Υ | 5 / Week | Bag length 28cm (11") Transparent | 0- |
| SOS-0039 | 2 Piece SenSura Mio Click Drainable Pouch Blue | SenSura Mio Click Maxi | Coloplast 11483 | Coupling/Flange 60mm (2 3/8") | Y | Maximum 10 | Easiclose™ Wide Outlet | |
| SOS-0040 | 2 Piece SenSura Mio Click Drainable Pouch Yellow | SenSura Mio Click Maxi | Coloplast 11487 | Coupling/Flange 70mm (2 3/4") | Υ | per patient admission | | - |
| SOS-0023 | 2 Piece Hollister Drainable Pouch Green | New Image Lock N Roll | Hollister 18132 | Coupling/Flange 44mm (1 3/4") | Y | of any combination of sizes | May only be ordered with the Hollister 2-Piece Barrier Bag length 30cm (12") | |
| SOS-0024 | IZ PIECE HOUISIEL DIAIDADIE POUCH REG | New Image Lock N Roll | Hollister 18133 | Coupling/Flange 57mm (2 1/4") | Y | | Transparent No filter | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | ΑΡΟ | Weekly Maximum | Notes | |
|-------------------------------|--|------------------|----------------------|--|----------|---------------------------|---|-----|
| One Piece | Urostomy System: Includes the cut-to-fit osto | my skin barrier | and urostomy pouc | ch as a single unit to collec | ct urine | with an integ | rated drain valve at the bottom to empty the bag. | |
| SOS-0041 | 1 Piece Urostomy Pouch Convex Light 10-23mm Cut-to-Fit | SenSura Mio | Coloplast 16835 | Stoma Size 10-23mm (3/8-7/8") | Y | | | |
| SOS-0042 | 1 Piece Urostomy Pouch Convex Light 15-33mm Cut-to-Fit | SenSura Mio | Coloplast 16836 | Stoma Size 15-33mm (5/8-1 1/4") | Y | 5 / Week | Light Convexity Bag length 26cm (10 1/4") Transparent | |
| SOS-0043 | 1 Piece Urostomy Pouch Convex Light 10-43mm Cut-to-Fit | SenSura Mio | Coloplast 16837 | Stoma Size 10-43mm (3/8-1 3/4") | Y | Maximum 10 per patient | | 2 |
| SOS-0044 | 1 Piece Urostomy Pouch Convex Deep 10-23mm Cut-to-Fit | SenSura Mio | Coloplast 16865 | Stoma Size 10-23mm (3/8-7/8") | Y | admission of any | | |
| SOS-0045 | ICul-IO-FIL | SenSura Mio | Coloplast 16866 | Stoma Size 15-33mm (5/8-1 1/4") | Y | combination of sizes | Deep Convexity Bag length 26cm (10 1/4") Transparent | 37 |
| SOS-0046 | 1 Piece Urostomy Pouch Convex Deep 10-43mm Cut-to-Fit | SenSura Mio | Coloplast 16867 | Stoma Size 10-43mm (3/8-1 3/4") | Y | | | |
| One Piece until it require | • | ble pouch with a | a larger opening for | ileostomies and colostom | nies. Ar | n integrated ro | Il closure or clamp is used to keep the pouch closed | |
| SOS-0047 | 1 Piece Drainable Pouch Flat 10-55mm Cut-to- Fit | SenSura Mio | Coloplast 10481 | Stoma Size 10-55mm (3/8-1 1/8") | Y | | | |
| SOS-0048 | 1 Piece Drainable Pouch Convex Light 10-23mm Cut-to-Fit | SenSura Mio | Coloplast 16737 | Stoma Size 10-23mm (3/8-7/8") | Y | | - 0 |)-) |
| SOS-0049 | 1 Piece Drainable Pouch Convex Light 15-33mm Cut-to-Fit | SenSura Mio | Coloplast 16738 | Stoma Size 15-33mm (5/8-1 1/4") | Y | 5 / Week | Bag length 28cm (11") | |
| SOS-0050 | 1 Piece Drainable Pouch Convex Light 10-43mm Cut-to-Fit | SenSura Mio | Coloplast 16739 | Stoma Size 10-43mm (3/8-1 3/4") | Y | Maximum 10 per patient | Wide Outlet | F |
| SOS-0051 | 1 Piece Drainable Pouch Convex Deep 10-23mm Cut-to-Fit | SenSura Mio | Coloplast 16755 | Stoma Size 10-23mm (3/8-7/8") | Y | admission of any | Transparent or Opaque (with inspection window) | |
| SOS-0052 | 1 Piece Drainable Pouch Convex Deep 15-33mm Cut-to-Fit | SenSura Mio | Coloplast 16756 | Stoma Size 15-33mm (5/8-1 1/4") | Y | combination of sizes | | |
| SOS-0053 | 1 Piece Drainable Pouch Convex Deep 10-43mm Cut-to-Fit | SenSura Mio | Coloplast 16767 | Stoma Size 10-43mm (3/8- 1 3/4") | Y | | | |
| SOS-0071 | 1 Piece High Output Pouch 10-100mm Cut-to-Fit | SenSura Mio | Coloplast 18666 | Baseplate hole size: 10-100mm - 725ml | Y | | Trasparent cut-to-fit 10-100mm | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|--|-------------|----------------------|--------------------------------|-----|--|--|
| | Piece System: Two-piece system with adhesive months to four years. The pediatric system is sire | | | | | | uring freedom of movement while still wearing a two-piece. For children e applied without pressure on the peristomal area. |
| SOS-0069 | Kids Flex Barrier Flat Cut-to-Fit | SenSura Mio | Colplast 18747 | Baseplate hole size: 0-35mm | Y | 5 / Week Maximum 10 per patient | Cut- to- fit 0-35 mm Designed for kids aged six months to four years, to be worn in combination with a compatible SenSura Mio Kids Drainable Pouch SSO-0070 |
| SOS-0070 | Kids Flex Drainable Pouch Midi | SenSura Mio | Coloplast 18724 | Coupling Size 35 mm | Y | admission of any combination of sizes | Opaque 35 mm Designed for kids aged six months to four years, to be worn in combination with a compatible SenSura Mio Kids barrier SSO-0069 |
| | | ·, | | | | | |
| SOS-0054 | Ostomy Support Belt for SenSura Mio Standard | Brava Belt | Coloplast 0423 | Each up to 100 cm (40") | Y | 1 per Patient | Standard Size |
| SOS-0055 | Ostomy Support Belt for SenSura Mio XXL | Brava Belt | Coloplast 0424 | Each up to 135 cm (61") | Y | Admission | XXL Size |
| SOS-0056 | Ostomy Strip Paste | Brava | Coloplast 2655 | Each 57g (2 oz) | Y | 5 / Patient Admission | Moldable. Fills gaps to create a tight seal to prevent leakage. |
| SOS-0057 | Ostomy Protective Sheet | Brava | Coloplast 3210 | Each 10 x10 cm 4 x 4" | Y | 3 | Absorbs moisture and protects peristomal skin when used under the barrier to minimize risk of maceration |
| SOS-0058 | Ostomy Paste (Tube) | Brava | Coloplast 12050 | Each Tube 60g (2 oz) | Y | 1 / Patient Admission | Alcohol-free paste fills folds and creases to maintain a tight fit between the stoma and the ostomy barrier to reduce leakage |
| SOS-0059 | Ostomy Powder | Brava | Coloplast 1907 | Each Bottle 25g (1 oz) | Y | 1 / Patient Admission | It absorbs moisture so the skin is kept dry to reduce irritation. Can be used around the stoma or under the barrier baseplate. |
| SOS-0066 | Mouldable Skin Barrier Ring | Brava | Coloplast 12045 | 18 mm diameter 4.2 mm thick | Y | 1 | A durable ring for reduced leakage. Creates a durable seal between the stoma and the base plate, thus reducing leakage. It is easy to use and can be molded to fit the patient's body shape. |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|---|-------------------|-----------------------|---|-----|-------------------|---|
| SOS-0065 | Ostomy Care Kit | Coloplast | Coloplast | Each Kit | Y | 1 | Maximum Quantity: 1 Kit/Patient Admission |
| One piece | e ostomy flat system | | | 1 | | | Coloplast |
| One piece | e ostomy convx light system | | | 1 | | | T |
| Scissors | | | | 1 | | | 6 |
| Mirror | | | | 1 | | | |
| Adhesive | Remover Wipe | | | 1 | | | |
| Lubricatir | ng Deordorant | | | 1 | | | Stonauto |
| Measurin | g Guide & Marker | | | 1 | | | |
| Coloplast | Care Information and Forms (Billingual) | | | 1 | | | Cotoplast |
| SOS-0060 | Fistula and Wound Pouch | Eakin | ConvaTec 839264 | For wounds up to 245 x 160 mm (9.7" x 6.3") | Y | 5 | Flexible pouches designed for skin protection and contain drainage from wounds, fistulas and difficult stoma sites. Horizontal cutting area. Suitable with remote drainage attachment, tap closure & window for wound access. |
| SOS-0061 | Skin Barrier Ring Small | Eakin Cohesive | ConvaTec 839002 | Each 50 mm (2") | Y | 1 | Can be moulded to a variety of shapes to help prevent leakage around the stoma by acting as filler |
| SOS-0062 | Skin Barrier Ring Large | Eakin Cohesive | ConvaTec 839001 | Each 100 mm (4") | Y | Ι | for uneven skin surfaces. |
| SOS-0063 | Antimicrobial Ostomy Ring | Hydrofera Blue | Hydrofera HBRF2650 | 60 mm (2.5") Diameter | Y | 3 | Not to be used longer than 6 months.Not to be use on pregnant or nursing patients. |
| SOS-0064 | Ceramide Skin Barrier Ring | Adapt CeraRing | Hollister 8805 | 48 mm (2") Diameter | Y | 1 | Infused with ceramide to help protect the peristomal skin's natural moisture barrier to prevent leakage and dryness. |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|---------------------|------------|----------------------|-----------------------|-----|--------------------------|---|
| SOS-0067 | Silver Ionic Powder | AG+ | Medline MSC9405 | Each 5 gram Bottle | Y | 1 / Patient Admission | Ideal for wounds that are difficult to dress; Manages bacterial burden; Powder can convert any other dressing to an antimicrobial; Contraindications: third- degree burns, patients with a known sensitivity to silver and as a surgical implant; Do not use topical antibiotics in conjunction with this powder |

URINARY & DRAINAGE SUPPLIES

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|--|------------------|-----------------------|-------------------------|------------|-------------------|--|
| and Commun | nity Care Support Servcies North East-funded supplies ma | ay be authorized | only while the patier | nt is in need of nursin | ng servi | ces. ODSP, s | al surgery, or who have a condition that makes urination difficult. Home social assistance, private insurance, DVA must be accessed as possible ad coatings are available for special order only as per the Exception Item |
| diameter was | used to facilitate passage of blood clots, the size can be | decreased by 2- | 4 Fr per catheter cha | ange until the optima | al size is | reached. Lik | e the smallest size diameter that permits good urine flow. If a larger weise, balloon sizes should be small, with 5 mL balloons, reserving the requent blockage may need catheter changed every 2-3 weeks or more |
| SUR-0001 | Foley 2 Way Hydrogel Coated Latex 12Fr 5 mL | BARD Canada | BARD 365712 | 12 Fr 5 mL Balloon | Y | 2/month | |
| SUR-0002 | Foley 2 Way Hydrogel Coated Latex 14Fr 5 mL | BARD Canada | BARD 365714 | 14 Fr 5 mL Balloon | Υ | 2/month | |
| SUR-0003 | Foley 2 Way Hydrogel Coated Latex 16Fr 5mL | BARD Canada | BARD 365716 | 16 Fr 5 mL Balloon | Y | 2/month | |
| SUR-0004 | Foley 2 Way Hydrogel Coated Latex 18Fr 5 mL | BARD Canada | BARD 365718 | 18 Fr 5 mL Balloon | Y | 2/month | |
| SUR-0014 | Foley 2 Way Hydrogel Coated Latex 14Fr 30mL | BARD Canada | BARD 366714 | 14 Fr 30 mL Balloon | Y | 2/month | |
| SUR-0015 | Foley 2 Way Hydrogel Coated Latex 16Fr 30 mL | BARD Canada | BARD 366716 | 16 Fr 30 mL Balloon | Y | 2/month | |
| SUR-0016 | Foley 2 Way Hydrogel Coated Latex 18Fr 30 mL | BARD Canada | BARD 366718 | 18 Fr 30 mL Balloon | Υ | 2/month | |
| SUR-0017 | Foley 2 Way Hydrogel Coated Latex 20Fr 30 mL | BARD Canada | BARD 366720 | 20 Fr 30 mL Balloon | Υ | 2/month | |
| SUR-0116 | Foley 2 Way Red Latex Tiemann Coude Tip 14Fr 5 mL | BARD Canada | BARD 0102L14 | 14 Fr 5 mL Balloon | Υ | 2/month | For difficult insertions, coude-tip Foley catheters feature a bent tip that maneuvers around obstructions leading to the bladder (eg. men with |
| SUR-0117 | Foley 2 Way Red Latex Tiemann Coude Tip 16Fr 5 mL | BARD Canada | BARD 0102L16 | 16 Fr 5 mL Balloon | Υ | 2/month | enlarged prostate) |
| SUR-0027 | Foley 2 Way 100% Silicone 12 Fr 10mL | Folysil | Coloplast AA6112 | 12 Fr 10 mL Balloon | Υ | 2/month | 100% Silicone; Sterile; Single-patient use; Non-latex |
| SUR-0028 | Foley 2 Way 100% Silicone 14Fr 10mL | Folysil | Coloplast AA6114 | 14 Fr 10 mL Balloon | Υ | 2/month | 100% Silicone; Sterile; Single-patient use; Non-latex |
| SUR-0037 | Foley 2 Way 100% Silicone 16Fr 15mL | Folysil | Coloplast AA6116 | 16 Fr 15 mL Balloon | Y | 2/month | 100% Silicone; Sterile; Single-patient use; Non-latex |
| SUR-0046 | Foley 2 Way 100% Silicone 18Fr 15mL | Folysil | Coloplast AA6118 | 18 Fr 15 mL Balloon | Υ | 2/month | 100% Silicone; Sterile; Single-patient use; Non-latex |
| SUR-0033 | Foley 2 Way 100% Silicone 18Fr 30mL | Folysil | Coloplast AA6C18 | 18 Fr 30 mL Balloon | Υ | 2/month | 100% Silicone; Sterile; Single-patient use; Non-latex |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|--|---------------------------|------------------------|----------------------------|------------|-------------------|---|
| SUR-0034 | Foley 2 Way 100% Silicone 20Fr 30mL | Folysil | Coloplast AA6C20 | 20 Fr 30 mL Balloon | Y | 2/month | 100% Silicone; Sterile; Single-patient use; Non-latex |
| SUR-0047 | Foley 2 Way 100% Silicone Coude Tip 14Fr 15mL | Folysil | Coloplast AA6314 | 14 Fr 15 mL Balloon | Y | 2/month | Tiemann (coude) tip for difficult insertions. Features a bent tip that maneuvers around obstructions leading to the bladder (eg. men with |
| SUR-0121 | Foley 2 Way 100% Silicone Coude Tip 16Fr 15mL | Folysil | Coloplast AA6316 | 16 Fr 15 mL Balloon | Y | 2/month | enlarged prostate). Latex free |
| INTERMITTE | | rile intermittent ca | theterization. All-ir | -one ready to use se | elf-lubica | ating catheter | and sterile bag set. The unique closed system helps to reduce the ris |
| SUR-0112 | Female Intermittent Catheter 10 Fr | SpeediCath Compact Set | Coloplast 28520 | 10 Fr x 9 cm | Y | 35 | |
| SUR-0113 | Female Intermittent Catheter 12 Fr | SpeediCath Compact Set | Coloplast 28522 | 12 Fr x 9 cm | Y | | Bag is directly attached to the catheter as a closed system Single patient use - Disposable. |
| SUR-0114 | Female Intermittent Catheter 14 Fr | SpeediCath Compact Set | Coloplast 28524 | 14 Fr x 9 cm | Υ | 35 | Available in male and female versions Maxiumum 5/day |
| SUR-0115 | Male Intermittent Catheter 12-18 Fr | SpeediCath Compact Set | Coloplast 28422 | 12-18 Fr 30 cm long | Υ | 35 | |
| SUR-0125 | Male External Condom Catheter Silicone 23mm | Freedom Clear | Coloplast 505711 | 23 mm | Y | 7 | Standard Adhesive; Non-sterile; Single Patient Use; Disposable |
| SUR-0126 | Male External Condom Catheter Silicone 28mm | Freedom Clear | Coloplast 505731 | 28 mm | Y | 7 | Standard Adhesive; Non-sterile; Single Patient Use; Disposable |
| SUR-0127 | Male External Condom Catheter Silicone 31mm | Freedom Clear | Coloplast 505741 | 31 mm | Y | 7 | Standard Adhesive; Non-sterile; Single Patient Use; Disposable |
| | RAINAGE BAGS: Urine drainage bags are used to c) are available by an approved Exception Item Requ | | | | use leg | bags are fo | r use with Foley Catheters. Drainage bags with slide valve (instea |
| | Leg Bag Kit 600 mL Flip-Flo Valve with Straps and 18" Tubing | Bard Canada | Dispoz-A-Bag 150180 | 19oz / 562ml 18" Tubing | Y | 2 | - For long-term catheter users, the bags should be replaced every 5-7 days |
| SUR-0086 | Leg Bag 562 mL Flip-Flo Valve with Straps | Bard Canada | Dispoz-A-Bag 1501MF | 19 oz / 562ml Each | Y | 2 | - Bags must be emptied when 1/3 to 1/2 full - Adjustable straps are included and attached to the leg bags. |
| SUR-0089 | Leg Bag 946mL Flip Flo Valve with Straps | Bard Canada | Dispoz-A-Bag 1501LF | 32 oz / 946ml Each | Y | 2 | - No longer best practice to change bag for overnight use unless the volume requires this to prevent leakage/overflow; |
| | Drainage Bag 2000mL with Anti-Reflux Device, Catheter tubing and Metal Drain Clamp | Bard Canada | Bard 154002 | 2000 mL Each | Y | 2 | - Larger sized 2000 ml bag often used for patients who are bedridden or end of life. |
| SUR-0094 | Urinary Leg Bag Extension Tubing 18" | Medline | Medline DYND12550 | 18 inches Each | Y | 2/month | 18" extension tubing can be cut to preferred length to meet patient needs |
| SOS-0013 | Nephrostomy Drainage Bag 600 mL | Argon | DBAG600 | 600 mL | Y | 1 | With adsjustable elastic strap for waist or leg. Provide 2 at admission then 1 / week until discharge. |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | ΑΡΟ | Weekly Maximum | Notes |
|---------------------|--|-------------------------|--------------------------------|---|---------|-------------------|---|
| CATHETER | STABILIZATION: Adhesive devices used to stabiliz | e indwelling urina | ary catheters, redu | cing catheter mover | nent a | nd minimizii | ng accidental dislodgments |
| SUR-0111 | Universal Foley Catheter Stabilization Device | StatLock | Bard FOL0102 | Each | Υ | 2/month | Maximum for nephrostomy care is up to 4/mth |
| SUR-0104 | Statlock Universal Plus Catheter Stabilization Device Large 12-14 FR | StatLock | Bard VUPD1214 | 12 - 14 fr | Y | 2/month | |
| SUR-0105 | StatLock Universal Plus Catheter Stabilization Device Medium 10 - 12 FR | StatLock | Bard VUPD1012 | 10 - 12 fr | Y | 2/month | Restricts catheter motion and fits virtually all catheter sizes. For biliary and nephrostomy drainage and peritoneal lavage catheters. |
| SUR-0106 | StatLock Universal Plus Catheter Stabilization Device Small 6 - 8.5 FR | StatLock | Bard VUPD68 | 6 - 8.5 fr | Y | 2/month | |
| CATHETER | ACCESSORIES: There are various accessories that | it can be used to | assist with either i | ndwelling catheter u | se or t | he process | of intermittent catheterisation. |
| SUR-0078 | Catheter Plug and Drain Tube Cover Set | Medline | Medline DYND12200 | Each | Y | 2/month | Fits securely into funnel end of most Foley catheters; Sterile; Cover protects connector on drainage bag |
| SUR-0110 | Syringe Catheter Tip 60 mL | BD | BD 309620 | 60 mL 2 oz tip | Y | 7 | Used for flushing/cleaning catheters, gastronomy tubes and other devices |
| SUR-0120 | Bulb Irrigation Syringe 60 mL | Medline | Medline DYND20125 | 60 mL | Y | 7 | Soft, pliable bulb for maximum control of drawn or expelled volume. Easy to grip and use with only one hand; Sterile |
| SUR-0091 | Sterile Bacteriostatic Lubricant Water Based | Healthcare Plus | Healthcare Plus 248-LUB-035 | 3.5 g pouch | Y | 2/month | For use only with Foley catheters. Not required for SpeediCath lubricated intermittent catherters. |
| SUR-0103 | Drainage Tube Universal 1 Way Stopcock | MED-Rx | MED-Rx 10-1210 | 16 inch | Y | 2/month | Maximum for nephrostomy care is up to 4/mth |
| SUR-0108 | Specimen container | Star Plex Scientific | Fisher Scientific 360993197 | 90 mL | Y | 2 | Sterile polypropylene container for collection, transport and storage of specimens, Single use, Disposable |
| | US BLADDER IRRIGATION (CBI): Approved for times the set of the set | e-limited treatm | ent for Continuous | Bladder Irrigation of | nly. U | sually requir | ed to discharge patient from hospital for end of life or hospice care. |
| SUR-0122 | 0.9% Sodium Chloride Irrigation Solution 3000 mL | Baxter | Baxter JB7127I | 3000 mL | N | by | Must be authorized/ordered by Care Coordinator SPECIAL ORDER ITEM: Availability of stock with the vendor needs to be confirmed prior to ordering. |
| SUR-0123 | Bladder Irrigation Tubing Set | Baxter | Baxter JC4040 | 0.195" Inner Dia. Length: 2.1m (81″ Long) | Y | 1 | Includes tubing, drip chamber, regulating clamp, and catheter adapter. Plug directly into 3-way catheter or us catheter adaptor unit when 2- way catheter is used. |
| SUR-0124 | Foley 3 Way Catheter for CBI Silicone 22 Fr | X-Flow | Coloplast AB6122 | 22 Fr 30-50 mL Balloon | Y | 1 | Patient usually discharged from hospital with required catheter. Couvalier tip appropriate for both male and female patients. Additional patient-specific sizes available with Exception Item Request approval |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes |
|---------------------|--|------------------|----------------------|----------------------|---------|-------------------|--|
| | ITS: Disposable, Single use prepared kits of commutional items may also be ordered as needed. | only used items. | Most items in eac | h kit can be ordered | indivio | lually if an e | ntire kit is not required if an individual CHRIS order code is listed |
| SUR-0079 | Foley Catheter Tray | MED-Rx | CHS 85-5050 | Each Tray | Y | 2 / Month | For catherter insertion 1/month |
| | CSR Wrap Blue | | | 1 | | | Order catheters, leg bags or drainage bags seperately |
| | Plastic Transfer Forcep | | | 1 | | | (may require 1 - SIV-0082 - 10mL Luer Lok syringe to deflate balloon for removal of catheter) |
| | Graduated Basin | | | 1 | | | If irrigation required add: |
| | Gloves Medium (Pair) Non-latex | | | 1 | | | SUR-0080 Irrigation Tray Kit with 60 mL Piston Syringe and SSO-0019 Normal Saline 500mL bottle |
| | Fenestrated Drape Paper | | | 1 | | | |
| | Towel Paper | | | 1 | | | |
| | Povidone lodine Swabstick (3 pack) | | | 2 | | | |
| | Lubricant Water Based | | | 1 | | | |
| | Foley Catheter Stabilization Device | | | 1 | | | |
| | Sterile Water 10 mL Pre-Filled Syringe | | | 1 | | | |
| SUR-0080 | Irrigation Tray Kit with 60 mL Piston Syringe | MED-Rx | CHS 85-5042 | Each Kit | Υ | 1 | Disposable, single use |
| | Irrigation Tray | | | 1 | | | |
| | 500 mL Graduated Container | | | 1 | | | |
| | 60 mL Piston Irrigation Syringe with Cap |] | | 1 | | | |
| | Alcohol Prep Pad | | | 1 | | | |
| | PVP Prep Pad | | | 1 | | | |
| | Underpad 17" x 18" | | | 1 | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | ΑΡΟ | Weekly Maximum | Notes |
|---------------------|---|---------------------------------|----------------------------|----------|-----|-------------------|---|
| SMK-0008 Pe | ercutaneous Nephrostomy Kit: | NE LHIN | Various - See each item | Each Kit | Y | 1 | Patient usually discharged from hospital with required catheter. Additional patient-specific sizes available as per this catalog or by |
| SDR-0005 | Dressing Tray Kit - Disposable Single Use | | | 2 | | | Exception Item Request approval |
| SDR-0055 | Drain Gauze Sponge Trach Sterile Woven 4"x4" | | | 2 | | | |
| SOT-0001 | Disposable Ear Loop Surgical Face Mask | | | 2 | | | |
| SOT-0021 | Sterile Nitrile Gloves Medium (Pair) | | | 2 | | | |
| SSO-0021 | Normal Saline 10 mL Pre-Filled Syringe | | | 1 | | | |
| SSO-0052 | Alcohol 70% Wipe | | | 14 | | | |
| SUR-0085 | Leg Bag Kit 600 mL Flip-Flo Valve with Straps and 18" Tubing | | | 1 | | | |
| SWC-0128 | Transparent Film Dressing 15x20cm | | | 2 | | | |
| SUR-0103 | Drainage Tube Universal 1 Way Stopcock | | | 2 | | | |
| SMK-0044 Pa | aracentesis Drainage Kit | | | | Y | 1 | To be ordered by and provided to Nurse Practitioners or Physicians only for abdomen or chest drainage procedures for eligible palliative |
| Not Applicable | Paracentesis Drainage Tray w/9cm Catheter | Smartap | CardioMed CM-0020T | 1 | | | patients to support end of life care. Home and Community Care Support Services North East will provide |
| SDR-0086 | White Surgical Tape 1" Single Use Roll | | | 1 | | | the medical supplies to the physician for eligible patients only. Funding or supplies for topical or local anesthetic for drainage procedures is not available and is to be provided by the physician or is |
| SOT-0001 | Disposable Ear Loop Surgical Face Mask | ble Ear Loop Surgical Face Mask | | | | | at the cost of the patient. |
| SCM-0019 | Blue Absorbent Pad 23"x24" | | | 2 | | | |
| SSO-0006 | Chlorhexidine 2% Alcohol 70% Swabstick | | | 2 | | | |
| SOT-0021 | Sterile Nitrile Gloves Medium (Pair) | | | 2 |] | | |

FEEDING, SUCTION & TRACHEOSTOMY SUPPLIES

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | ΑΡΟ | Weekly Maximum | Notes | | | | |
|---------------------|--|-----------------------|---------------------------------|--------------------------------|--------|------------------------|---|--|--|--|--|
| feeding supp | itial supply authorization for four weeks only for short-term treatment and to allow patients to coordinate insurance or ADP applications. Supplies will only continue to be provided as per approved policy. Funder eding supplies shall only be provided for Kangaroo ePump or Joey Pump and may not be compatible with other pumps. Home and Community Care Support Services North East is not obligated to provide applies beyond defined periods while patients wait for third-party funding. | | | | | | | | | | |
| ENTERAL | FEEDING | | | | | | | | | | |
| ETM-1002 | auto flush functionality NOTE: Equipment is available for rental as per approved policies from Medigas (Praixair). | | | | | | | | | | |
| ETM-1005 | IV pole - portable | Please review | to the Medigas Med | dical Equipment Ca | talogs | for more detai | led information before ordering. | | | | |
| SFS-0005 | Enteral Feeding Gravity Bag | Kangaroo | Cardinal Health Z8884702500 | 1000 ml | Y | 2 | To be used in the absence of a feeding pump. Non-sterile | | | | |
| SFS-0006 | Enteral Feeding Pump Set 500 ml | Kangaroo Joey | Cardinal Health Z762055 | 500 ml | Y | 2 | | | | | |
| SFS-0031 | Enteral Feeding Pump Set 1000 ml | Kangaroo Joey | Cardinal Health Z763656 | 1000 ml | Y | 2 | | | | | |
| SFS-0035 | Enteral Feeding Pump Set with Flush Set 1000ml | Kangaroo Joey | Cardinal Health Z763662 | 1000 ml | Y | 2 | For continuous overnight feeds ONLY | | | | |
| SFS-0045 | Feeding Tubing Set with Safety Screw Spike | Kangaroo Joey | Cardinal Health Z765559 | 1 Tube | Y | 7 | By request from Physician or Registed Dietician only for patients with a clinical need for a closed system. Tubing only - no bag. | | | | |
| SFS-0010 | Gastrostomy Tube with Y Port 16 Fr | Kangaroo | Cardinal Health Z8884720163E | 16 Fr x 9" 20 cc | Y | 1/Patient Admission | | | | | |
| SFS-0036 | Gastrostomy Tube with Y Port 20 Fr | Kangaroo | Cardinal Health Z8884720205E | 20 Fr x 9" 20 cc | Y | 1/Patient Admission | | | | | |
| SFS-0012 | Nasogastric Feeding Tube 5 Fr (Infant) | Kangaroo | Cardinal Health Z461412E | 5 Fr x 20" 1.7 mm x 51 cm | Y | 1/month | Remains in situ for 30 days Purple Polyurethane Sterile | | | | |
| SFS-0037 | Nasogastric Feeding Tube PVC 5 Fr (Infant) | Kangaroo | Cardinal Health Z460802E | 5 Fr x 16" 1.7 mm x 41 cm | Y | 2/week | Used only for patients who pull them out frequently PVC Sterile | | | | |
| SFS-0038 | Nasogastric Feeding Tube 6.5 Fr (Infant) | Kangaroo | Cardinal Health Z461438E | 6.5 Fr x 36" 2.2 mm x 91 cm | Y | 1/month | Remains in situ for 30 days Purple Polyurethane Sterile | | | | |
| SFS-0015 | Nasogastric Feeding Tube 12 Fr | Kangaroo Entriflex | Cardinal Health Z8884721252E | 12 Fr x 43" 4 mm x 109 cm | Y | 1/month | With stylet and rigid outlet port | | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | ΑΡΟ | Weekly Maximum | Notes |
|---------------------|---|----------------------|---------------------------------|-----------------------------|-----|------------------------|--|
| SFS-0040 | Adapter Y-Port-G-tube/P.E.G. 16 Fr | Kangaroo | Cardinal Health Z8884751622E | Fits 16 Fr or smaller | Y | 1/Patient Admission | Non-Sterile. For use with the Kangaroo P.E.G. Kits and can be used as the replacement Y-port on all Kangaroo gastrostomy tubes. 1 per patient |
| SFS-0041 | Adapter Y-Port-G-tube/P.E.G. 20 Fr | Kangaroo | Cardinal Health Z8884752042E | Fits all 18 Fr or bigger | Y | 1/Patient Admission | admission. Replace only if port is clogged or broken when patient still eligible for supplies. |
| SFS-0044 | Feeding Tube Y-Site Extension Set | Kangaroo | Cardinal Health Z8884705008 | 12 inch | Υ | 1 | Not compatible with low-profile tubes like MIC-KEY or Nutriport |
| SFS-0042 | Purple Enteral Oral Syringe 3 mL | Monoject | Cardinal Health Z403S | 3 mL | Y | 3 | Standard Tip; Non sterile; Can be reused if washed with same patient For Neonatal and Pediatric Feeding Systems |
| SFS-0043 | Purple Enteral Oral Syringe 10 mL | Monoject | Cardinal Health Z410S | 10 mL | Y | 3 | Standard Tip; Non sterile; Can be reused if washed with same patient For Neonatal and Pediatric Feeding Systems |
| SES-0047 | Purple Enteral Oral Syringe 35 mL Sterile with Tip Cap | Monoject | Cardinal Health Z435SG | 35 mL | Y | 3 | Standard Tip; Sterile; For Neonatal and Pediatric Feeding Systems Single Patient Use with Nasogastric Feeding Tube |
| SFS-0048 | 3 mL Enteral Syringe with ENFit Connection | Monoject ENFit | Cardinal Health Z403SE | 3 mL | Υ | 3 | Can be reused if washed with same patient ENFit® Connection |
| SFS-0049 | 12 mL Enteral Syringe with ENFit Connection | Monoject ENFit | Cardinal Health Z412SE | 12 mL | Υ | 3 | Can be reused if washed with same patient ENFit® Connection |
| SFS-0050 | 35 mL Enteral Syringe with ENFit Connection | Monoject ENFit | Cardinal Health Z435SE | 35 mL | Υ | 3 | Can be reused if washed with same patient ENFit® Connection |
| SFS-0051 | 60 mL Enteral Syringe with ENFit Connection | Monoject ENFit | Cardinal Health Z460SE | 60 mL | Υ | 3 | Can be reused if washed with same patient ENFit® Connection |
| SFS-0046 | Low Profile Tube Feed Extension Set | MIC-KEY Secur-Lok | Kimberly Clark 0121-12 | 12 inch | Y | 1 | Primarily used with pediatric patients. SECUR-LOK Right-Angle Connector; 2 Port "Y" connector and Tubing clamp for use with low- profile MIC-KEY tubes. Not compatible with Nutriport low-profile tubes. |
| SFS-0039 | Nasogastric Tube Stabilization Device (Pediatric) | Statlok | Bard NG0100 | Pediatric | Y | 1 | Order for babies and small children For nasogastric tubes and feeding tubes |
| SFS-0019 | Nasogastric Tube Stabilization Device (Adult) | Statlok | Bard NG0120 | Adult | Y | 1 | For nasogastric tubes and feeding tubes |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | ΑΡΟ | Weekly Maximum | Notes | | | |
|---------------------|---|-----------------|----------------------|--------------------|--------|-------------------|--|--|--|--|
| SUCTION | SUCTION CATHETERS, TUBING AND ACCESSORIES | | | | | | | | | |
| ETM-1006 | Suction machine (includes nebulizer bottle, 2 trac masks, and tubing) NOTE: Equipment is available for rental as per approved policies from Medigas (Praixair). | | | | | | | | | |
| ETM-1015 | Gastric Drainage Aspirator | Please review t | o the Medigas Med | lical Equipment Ca | talogs | for more detai | iled information before ordering. | | | |
| SOT-0090 | Disposable Suction Tubing Kit for ETM-1015 | Gomco | Allied 01-90-2000 | 18" and 72" | Y | 7 | Includes two different length tubes per kit | | | |
| SOT-0091 | Disposable Filter (Box of 3) for ETM-1015 | Gomco | Allied 01-90-3100 | Box of 3 | Υ | 1 | Filter only needs replacing if it gets wet | | | |
| SOT-0092 | Disposable Collection Bottle for ETM-1015 | Gomco | Allied 01-90-3695 | 1100 mL | Y | 7 | Plastic; Single Patient Use | | | |
| SOT-0101 | Connector for Suction Tubing for ETM-1015- ONLY | Medline | DYND50510 | Each | Y | 7 | 5–in-1 straight plastic connector to connect tubing to ETM-1015 Gastric Drainage Aspirator | | | |
| SOT-0072 | Universal Suction Tubing with Connectors | Ornex | Medline OR66A | 1/4" X 72" | Y | 1 | Straight connectors feature ribbing which provides a firm connection to suction devices of all kinds. Designed to handle 28 Hg pressure. | | | |
| SOT-0053 | Suction Catheter 8 FR | Med-Rx | CHS 50-1008 | 8 Fr x 15" | Y | 7 | | | | |
| SOT-0054 | Suction Catheter 12Fr | Med-Rx | CHS 50-1012 | 12 Fr x 21" | Υ | 7 | Angled tip; Open; With control valve; Without tubing; Straight, Non- latex; Sterile | | | |
| SOT-0055 | Suction Catheter 14Fr | Med-Rx | CHS 50-1014 | 14 Fr x 21" | Y | 7 | | | | |
| SOT-0057 | Yankauer Suction Bulb Tip | Medline | DYND50130 | Each | Y | 1 | | | | |
| TRACHEC | STOMY PRODUCTS | | | | | | · | | | |
| | tional products are available in the TRACHEOT ated inner cannulas (4IC65; 5IC70; 6IC75; 7IC8 | | | · · · · · | - | • • • • | Compression & Wraps section; | | | |
| SOT-0032 | Tracheostomy Care Kit | Med-Rx | 85-4037 | Each | Y | 7 | Includes Disposable Tray, 4x4 Gauze, Plastic transfer Forcep, Pipe Cleaners, Cotton Tipped Applicators, Trach Sponges, Trach Brush, Twill Tape | | | |
| SOT-0034 | Tracheostomy Ties One Size | Shiley | Medtronic TTH | Each | Y | 1 | Fits both adults and pediatric patients with neck size up to 18" | | | |

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | ΑΡΟ | Weekly Maximum | Notes | |
|---------------------|--|------------------|----------------------|---|--------|-------------------|---|--|
| SOT-0014 | Trach Inner Cannula Disposable Beveled Tip for 6.5mm Tube Green | Shiley | Medtronic 4IC65 | I.D. x Length 5.5mm x 62mm | Y | 21 | For use with Tracheostomy Tubes Green 4CN65H or 4UN65H BPA & DEHP Free; Non-Latex; Non-sterile; Single-patient use | |
| SOT-0102 | Trach Inner Cannula Disposable Beveled Tip for 7.0mm Tube Orange | Shiley Flexible | Medtronic 5IC70 | I.D. x Length 6.0mm x 68mm | Υ | | For use with Tracheostomy Tube Orange 5CN70H or 5UN70H BPA & DEHP Free; Non-Latex; Non-sterile; Single-patient use | |
| SOT-0015 | Trach Inner Cannula Disposable Beveled Tip for 7.5mm Tube Light Blue | Shiley Flexible | Medtronic 6IC75 | I.D. x Length 6.5mm x 74mm | Υ | 21 | For use with Tracheostomy Tube Light Blue 6CN75H or 6UN75H BPA & DEHP Free; Non-Latex; Non-sterile; Single-patient use | |
| SOT-0103 | Trach Inner Cannula Disposable Beveled Tip for 8.0mm Tube Purple | Shiley Flexible | Medtronic 7IC80 | I.D. x Length 7.0mm x 77mm | Υ | 21 | For use with Tracheostomy Tube Purple 7CN80H or 7UN80H BPA & DEHP Free; Non-Latex; Non-sterile; Single-patient use | |
| SOT-0018 | Trach Inner Cannula Disposable Beveled Tip for 8.5mm Tube Navy Blue | Shiley Flexible | Medtronic 8IC85 | I.D. x Length 7.5mm x 79mm | Υ | 21 | For use with Tracheostomy Tube Navy 8CN85H or 8UN85H BPA & DEHP Free; Non-Latex; Non-sterile; Single-patient use | |
| SOT-0080 | Fenestrated Inner Cannula Disposable Size 4 | Shiley Legacy | Medtronic 4DICFEN | 4.9mm x 62mm | Υ | 21 | For use with Tracheostomy Tubes 4DFEN, 4DCFN | |
| SOT-0081 | Fenestrated Inner Cannula Disposable Size 6 | Shiley Legacy | Medtronic 6DICFEN | 6.0mm x 62mm | Y | 21 | For use with Tracheostomy Tubes 6DFEN, 6DCFN | |
| SOT-0082 | Fenestrated Inner Cannula Disposable Size 8 | Shiley Legacy | Medtronic 8DICFEN | 7.5mm x 79mm | Y | 21 | For use with Tracheostomy Tubes 8DFEN, 8DCFN | |
| The followin | g Tracheostomy tubes are to be ordered only v | when suctioning | g is required for F | ENESTRATED tra | cheost | tomy tubes i | n situ | |
| | Inner Cannula Disposable Size 4 for fenestrated tracheostomy tubes | Shiley Legacy | Medtronic 4DIC | 4.9mm x 62mm | Y | 21 | For use with fenestrated Tracheostomy Tubes 4DICFEN forsuctioning only | |
| $SOU_0/8$ | Inner Cannula Disposable Size 6 for fenestrated tracheostomy tubes | Shiley Legacy | Medtronic 6DIC | 6.0mm x 62mm | Υ | 21 | For use with fenestrated Tracheostomy tubes 6DICFEN for suctioning only | |
| ISOT_0079 | Inner Cannula Disposable Size 8 for fenestrated tracheostomy tubes | Shiley Legacy | Medtronic 8DIC | 7.5mm x 79mm | Y | 21 | For use with fenestrated Tracheostomy tubes 8DICFEN forsuctioning only | |
| OXYGEN I | MONITORING | | | | | | | |
| EPY-1008 | Oximetry Monitor | | | rental as per appro cal Equipment Cata | | | digas (Praixair). d information before ordering. | |
| SOT-0099 | Infant Toe Sensor for EPY-1008 O2 Monitor | Masimo Canada | Masimo 2328 | Each | Y | 1 | Single Patient Use. Includes one adhesive tape. | |
| SOT-0100 | Replacement Tapes for Infant Toe Sensor | Masimo Canada | Masimo 2307 | Each | Υ | 5 | Single Patient Use, Replacement adhesive tape for SOT-0099 | |

PROTECTIVE WEAR

| CHRIS Order Code | Description | Brand Name | Manufacturer Code | Size | APO | Weekly Maximum | Notes | |
|---|--|---|--|----------------------------------|--------------------|---------------------------------|--|--|
| hazards as pe provided PPE requires prote | IVE WEAR: Personal protective equipment (PPE) is a er policy. Service Providers are expected to provide the E if handling hazarous chemicals/drugs for approved pati ective wear be worn to minimize risk to patient as approve ents require a sterile technique as per approved policies | ir own personal pr ents. Non-Sterile ved by the Care Co | otective equipment PPE may be author pordinator. Clean P | to their own s ized for contr | taff to acted : | protect from i Service Provi | nfectious diseases. Service Providers may be ders <u>only</u> when the on-going patient-specific treatmer | |
| SOT-0004 | Disposable Isolation Gown Regular/Large | Medline | Medline CRI4000 | Each | Ν | 1 per Visit | Single use disposable; For use for protection from chemotherapy or other cytotoxic medications in liquid | |
| SOT-0098 | Disposable Chemo Gown Regular/Large | Medline | Medline NON27457 | Each | Ν | 1 per Visit | form (Intravenous, subcutaneous, intramuscular, liquoral, etc.). NOTE: This face mask is not appropriate | |
| SOT-0001 | Disposable Ear Loop Surgical Face Mask | Primaguard 120 | priMED Medical PG4-1092 | Each | Ν | 1 per Visit | the cytotoxic medication is aerosolized - Use a fit- tested N95 mask provided by the nurse's employer. | |
| SOT-0019 | Eye Protection Safety Glasses with Side Shields | Medical Mart | 320-NON24777V | Each | N | 1 | For procedures and patient-care activities likely to generate splashes or sprays of blood or body fluids or cytotoxic medications in liquid form Non-latex, powder free, nitrile examination gloves with 11 in. cuff. Approved for use with cytotoxic medications. Nursing providers are required to provide their staff | |
| SOT-0020 | Sterile Nitrile Gloves Small (Pair) | Flexam | Cardinal Health N8830 | Each Pair | Y | 3 | | |
| SOT-0021 | Sterile Nitrile Gloves Medium (Pair) | Flexam | Cardinal Health N8831 | Each Pair | Y | | | |
| SOT-0022 | Sterile Nitrile Gloves Large (Pair) | Flexam | Cardinal Health N8832 | Each Pair | Y | | with non-sterile gloves. Sterile gloves will <u>only</u> to be provided to nursing service provider personnel wher | |
| SOT-0023 | Sterile Nitrile Gloves Extra Large (Pair) | Flexam | Cardinal Health N8833 | Each Pair | Υ | 3 | sterile technique is required for patient care as approved by the Care Coordinator | |
| SOT-0012 | Sterile Nitrile Gloves Medium (Box) | Flexam | Cardinal Health N8831 | Box of 40 Pairs | N | 1 | Non-latex, powder free, nitrile examination gloves Bulk orders available Clinics and Residental Hospic only. Approved for use with cytotoxic medications. | |
| SOT-0093 | Non-Sterile Nitrile Gloves Small (Pair) | Esteem XP | Cardinal Health N8851XP | Each Pair | N | 7 | Non-latex, powder free, nitrile examination gloves. | |
| SOT-0094 | Non-Sterile Nitrile Gloves Medium (Pair) | Esteem XP | Cardinal Health N8852XP | Each Pair | Ν | 7 | ASTM F 739 Chemicals tested, but not for use with chemotherapy. Nursing providers are required to provide their staff with non-sterile gloves. | |
| SOT-0096 | Non-Sterile Nitrile Gloves Large (Pair) | Esteem XP | Cardinal Health N8853XP | Each Pair | N | 7 | | |
| SOT-0095 | Non-Sterile Nitrile Gloves Medium (Box) | Esteem XP | Cardinal Health N8852XP | Box of 50 | Ν | | Bulk orders only available for Home and Community Care North East clinical staff or Clinics or Hospices | |

SERVICES DE SOUTIEN À DOMICILE ET EN MILIEU COMMUNAUTAIRE Nord-Est

Regional Approved Medical Supplies List

Appendix A Depot/Drop-off Locations

Revised July 2022



*Note: Vendor will contact client when supplies are ready for pickup at depot or to coordinate private-pay patient home delivery requests

ALGOMA SUB-REGION (Questions or issues call: 705-949-1808 / 1-800-668-7705 Ext: 3269)

| Service Area | Depot/Drop-Off Location | Pick-Up Location | Service Type | Key Notes |
|-----------------|--|--|---|---|
| URBAN | Soo Medical | 656 Great Northern Road, Sault Ste. Marie, ON P6B 4Z9 705-574-0228 Fax: 705-574-0231 | Soo Medical (Robinson's) (IV) Equipment & Supplies, Medical Supplies | Patient will receive a phone call when supplies are ready |
| URBAN | The Medicine Shoppe | 316 Second Line West Sault Ste. Marie, ON P6C 2J5 705-253-0720/877-374-6938 Fax: 705-256-7656 | Medicine Shoppe's Infusion (IV) Equipment & Supplies | for pickup. |
| URBAN | ParaMed Clinic | 477 Queen Street East Suite 100 Sault Ste. Marie, ON, P6A 1Z5 705-759-0110 | Soo Medical (Robinson's) Medical Supplies | For ParaMed Clinic Patients only |
| URBAN | Algoma Residential Hospice Care (ARCH) | 229 Fourth Line W Sault Ste., Marie, ON P6A 0B5 705-942-1556 | Soo Medical (Robinson's) Medical Supplies & Medicine Shoppe's Infusion (IV) Equipment & Supplies | For Hospice Residents only |
| RURAL | Matthews Memorial Hospital | 1180 Richards Street, Richards Landing, ON POR 1J0 705-246-2570 | Soo Medical (Robinson's) Medical Supplies | |
| RURAL | Main Street Pharmacy | 156 Main Street, Thessalon, ON POR 1L0 705 842 2322 | Soo Medical (Robinson's) Medical Supplies Medicine Shoppe's Infusion (IV) Equipment & Supplies | |
| RURAL | Pit Stop Gas Bar & Confectionary | 3245 Hwy 17 East, Echo Bay, ON POS 1C0 705-248-3333 | Soo Medical (Robinson's) Medical Supplies | |
| RURAL | Goulais River Country Store | 2771 Hwy 17 North Goulais River, ON P0S 1E0 705-649-3811 | Soo Medical (Robinson's) Medical Supplies | |
| RURAL | Elliot Lake PharmaChoice | 267 Hillside Drive South Elliot Lake, ON P5A 1N7 705-848-0004 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | Home Hardware | 275 Imperial St Massey, ON POP 1P0 705-865-2135 | Robinson's Medical Supplies | |
| RURAL | Fenlon's Pharmacy | 88 Broadway Ave. Wawa, ON POS 1K0 705-856-2800 | Soo Medical (Robinson's) Medical Supplies | |
| RURAL | Patient Service Treatment Location | Patients living in the defined Hornepayne area are pre- authorized for home delivery as there is no available depot. | Soo Medical (Robinson's) Medical Supplies Medicine Shoppe's Infusion (IV) Equipment & Supplies | |
| RURAL | Patient Service Treatment Location | Patients living in the defined Blind River area are pre-authorized for home delivery because there is no available depot. | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |

Deliveries during extended hours: Only available for New Patients or Patient Resumes/ Hospital Discharges/End of Life

Private-Pay Patient Home Deliveries: Robinson's to coordinate for Patients when requested

*Note: Vendor will contact client when supplies are ready for pickup at depot or to coordinate private-pay patient home delivery requests

NIPISSING-TEMISKAMING SUB-REGION (Questions or issues call: Robinson's Pharmasave 705-495 4515 - Option 3. *NOTE:* For Sturgeon Falls medical supplies only call 705 566 5551 - Option 4)

| Service Area | Depot/Drop-Off Location | Pick-Up Location | Service Type | Key Notes |
|-----------------|--|--|--|--|
| URBAN | VON Ambulatory Nursing Clinic | 150 First Ave.West North Bay, ON P1N 3B9 705-472-8050 Fax: 705-495-6691 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | For VON Clinic Patients only |
| URBAN | Robinson's Pharmasave 705-495-4515 | 2547 Trout Lake Road North Bay, ON P1B 7S8 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| URBAN | CarePartners Ambulatory Nursing Clinic | 195 First Ave. West, North Bay, ON P1B 3B8 1-866-886-4172 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | For CarePartners Clinic Patients only |
| URBAN | Nipissing Serenity Hospice | 799 John Street, North Bay, ON P1B 8T2 705-499-6295 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | For Hospice Residents Only |
| RURAL | Pharmacie Nipissing Pharmacy, PharmaChoice | 26 Principale Street East Verner, ON P0H 2M0 705-581-2013 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | Mattawa IDA | 300 Main Street Mattawa, ON P0H 1V0 705-744-5981 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | South River Pharmacy | 311 Hwy 124 South River, ON P8A 1X0 705-386-0636 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | Loblaws Pharmacy (No Frills) | 12035 Hwy #17 Sturgeon Falls, ON P2B 2S6 705-753-5850 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | Clement's IDA | 497 Main St Powassan, ON P0H 1Z0 705-724-3523 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | Argyle Nursing Station | 11851 Hwy 522 Port Loring, ON P0H 1Y0 705-757-1717 | Robinson's Medical Supplies | |
| RURAL | Arnica Medical Pharmacy | 162 Huston St. Burks Falls, ON POA 1C0 705-382-2111 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | Sundridge Guardian Pharmacy | 107 Main St. Sunridge, ON P0A 1Z0 705-384-5335 | Robinson's Medical Supplies and Infusion (IV) Equipment & Supplies | |

Private-Pay Patient Deliveries: Robinson's to coordinate for Patients when requested

*Note: Vendor will contact client when supplies are ready for pickup at depot or to coordinate private-pay patient home delivery requests

NIPISSING-TEMISKAMING SUB-REGION (Questions or issues call Kirkland Lake office at 1-888-533-2222 Ext: 3269 or call Robinson's Pharmasave 705-495-4515 (Option 3)

| Servic e Area | Depot/Drop-Off Location | Pick-Up Location | Service Type | Key Notes |
|------------------|-----------------------------------|--|--|--|
| RURAL | Pharmasave BDR Drug Mart | 15 Station Rd. N Kirkland Lake, ON P3N 3J1 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | J&S Variety Store | 91 Pine Street Elk Lake, ON P0J 1G0 705-678-2246 | Robinson's Medical Supplies | |
| RURAL | Elk Lake Nursing Station | 71 Front Street Elk Lake, ON P0J 1G0 | Robinsons Infusion (IV) Equipment & Supplies | |
| RURAL | Larder Lake Co-Op | 25 Godfrey Street Larder Lake, ON P0K 1L0 705-643-2327 | Robinson's Medical Supplies | |
| RURAL | Larder Lake Pharmasave | 15 Godfrey Street Larder, Lake, ON P0K 1L0 | Robinsons Infusion (IV) Equipment & Supplies | |
| RURAL | Dogstar Service | 418 Moyneur Ave Matachewan, ON P0K 1M0 705 565 2262 | Robinson's Medical Supplies | |
| RURAL | Matachewan Nursing Station | 1 Moyneur Ave Matachewan, ON P0K 1M0 705-565-2351 | Robinson's Medical Supplies | |
| RURAL | CarePartners Nursing Clinic | 883361 Highway 65 N New Liskeard, ON P0J 1P0 1-866-886-4172 ext:5351 | Robinson's Medical Supplies and Infusion (IV) Equipment & Supplies | For CarePartner's Clinic Patients only |
| RURAL | Loblaw Pharmacy | 55 Scott Street New Liskeard, ON P0J 1P0 | Robinson's Medical Supplies and Infusion (IV) Equipment & Supplies | |
| RURAL | Latchford Town Office | 10 Main St. Latchford, ON P0J 1N0 705-676-2416 | Robinson's Medical Supplies | |
| RURAL | Temagami Family Health Team | 17 O'Connor Dr. Temagami, ON POH 2H0 705-569-3244 | Robinson's Medical Supplies and Infusion (IV) Equipment & Supplies | |
| RURAL | Marshall's Pharamsave | 43 Third Street Englehart, ON P0J 1H0 705-544-7878 | Robinson's Medical Supplies and Infusion (IV) Equipment & Supplies | |
| RURAL | Cote's Variety Store | 16 Tenth St. E Earlton, ON P0J 1E0 705-563-2022 | Robinson's Medical Supplies | |
| RURAL | Bayshore Specialty Rx | 70 Esna Park Drive Unit 11 Markham, ON L3R 6E7 905-474-0822/888-313-6988 | Bayshore Rx IV Supplies & Equipment | Next business day delivery from Markham |

Deliveries during extended hours: Only available for New Patients or Patient Resumes/ Hospital Discharges/End of Life Private-Pay Patient Deliveries: Robinson's to coordinate for Patients when requested

*Note: Vendor will contact client when supplies are ready for pickup at depot or to coordinate private-pay patient home delivery requests

COCHRANE SUB-REGION (Questions or issues call: 1-888-668-2222)

| Service Area | Depot/Drop-Off Location | Pick-Up Location | Service Type | Key Notes | |
|-----------------|---|---|--|--|--|
| URBAN | Bayshore Home Health Clinic | 202-119 Pine Street S Timmins, ON P4N 2K3 705-268-6088/866-227-3661 Fax: 705-268-9085 | Home and Community Care Support Services Medical Supplies & Bayshore Rx IV Supplies & Equipment | | |
| URBAN/ RURAL | Bayshore Specialty Rx | 70 Esna Park Drive Unit 11 Markham, ON L3R 6E7 905-474-0822/888-313-6988 | Bayshore Rx IV Supplies & Equipment | Next business day delivery from Markham | |
| RURAL | Iroquois Falls | Patients living in the defined Iroquois Falls, Matheson area are pre-authorized for home | Home and Community Care Support Services Medical Supplies Bayshore Rx IV Supplies & Equipment delivered to home | Stock Replenishment Only at VON Clinic Anson General Hospital 58 Anson Drive Iroquois Falls, ON P0K 1E0 | |
| RURAL | Matheson | delivery as there is no available depot. | (Bayshore Rx BodyGuard IV pump available at Iroquois Falls and Matheson clinic) | For VON Clinic Patients <u>only</u> Bingham Hospital 507 Eight Avenue (Basement) Matheson, ON POK 1N0 705-273-4500 | |
| RURAL | Smooth Rock Falls Pharmacy | 133 5th Street Smooth Rock Falls, ON POL 2B0 705-338-4470 | Home and Community Care Support Services Medical Supplies | | |
| RURAL | Cochrane | | Home and Community Care Support Services Medical Supplies Bayshore Rx IV Supplies & | Stock Replenishment Only at VON Clinic P.O Box 1750, 22E 5 th Street Cochrane, ON POL 1C0 877-513-2626 ext 7300 | |
| RURAL | Hearst | Patients living in the defined Cochrane, Hearst and Kapuskasing areas are pre- authorized for home delivery as there is no available depot. | Equipment delivered to home (Bayshore Rx BodyGuard IV pump available at Cochrane and Hearst Clinic) | Stock Replenishment Only VON Clinic PO Box 547 1403B Edward Street Hearst ON POL 1N0 705-362-6633 ext 7280 | |
| RURAL | Kapuskasing | | | Stock Replenishment Only VON Clinic 101 Progress Cr. Kapuskasing, ON P5N 3H5 | |
| RURAL | Missinaibi-Mattice Clinic | 500 Hwy 11 East Mattice, ON P0L 1T0 705-364-2161 | Home and Community Care Support Services Medical Supplies | | |
| RURAL | Home and Community Care Support Services Moosonee Office | PO Box 57, 48 Revillion St. Moosonee ON, POL 1Y0 705-336-2233 Fax: 705-336-2919 | Home and Community Care Support Services Medical Supplies | All Bayshore Rx IV Supplies | |
| RURAL | Cavalier Motel | 98 Monk Street Chapleau, ON P0M 1K0 705-864-0841 | Home and Community Care Support Services Medical Supplies | & Equipment is delivered to home | |
| RURAL | Jr Corner Store Gogama | 41 Poupore St, Gogama, ON, POM 1W0 705-894-2540 Attention: Gogama Nursing Clinic | Home and Community Care Support Services Medical Supplies | | |

Private-Pay Patient Deliveries for NE Home and Community Care Support Services Supplies: Call NE Home and Community Care Support Services for availability.

Note: All IV orders from Bayshore Rx are delivered directly to the patient's treatment location

*Note: Vendor will contact client when supplies are ready for pickup at depot or to coordinate private-pay patient home delivery requests

SUDBURY-MANITOULIN-PARRY SOUND SUB-REGION

(Questions or issues call: Robinson's Pharmasave (705-566-5551 Option 4)

| Service Area | Depot/Drop-Off Location | Pick-Up Location | Service Type | Key Notes |
|-----------------|--|---|---|---|
| URBAN | Health Care Pharmacy Pharmasave | 1276 Lasalle Blvd. Sudbury, ON P3A 1Y8 (705) 566-5592 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| URBAN | Bayshore Ambulatory Nursing Clinic | 2120 Regent Street, Unit 8 Sudbury, ON P3E 5S8 705-523-6668 | Robinson's Medical Supplies and IV Equipment & IV Supplies | For Bayshore Clinic Patients only |
| URBAN | Corner Clinic | 344 Elgin, Sudbury, ON P3E 3N9 705-673-3721 | Robinson's Medical Supplies | |
| URBAN | Daffodil Lodge | 41 Ramsey Lake Rd. Sudbury, ON P3E 5J1 | Robinson's Medical Supplies | For patients residing at Daffodil Lodge only |
| URBAN | ICAN (Haig) | 582 Haig Street Sudbury, ON P3C 1C4 | Robinson's Medical Supplies | For ICAN patients only |
| URBAN | ICAN (Mont Adam) | 117 Mont Adam Sudbury, ON P3B 1E1 | Robinson's Medical Supplies | For ICAN patients only |
| URBAN | Sudbury Maison McCulloch Hospice | 1028 South Bay Sudbury, ON P3E 6J7 705-674-9252 | Robinson's Medical Supplies and IV Equipment & IV Supplies | For Hospice Residents only |
| URBAN | Health Sciences North | 41 Ramsey Lake Rd Sudbury, ON P3E 5J1 | Robinson's Medical Supplies and IV Equipment & IV Supplies | Not for Patient pick-up |
| URBAN | Chelmsford Guardian Pharmacy | 46 Main St. Chelmsford, ON POM 1L0 705-855-2371 | Robinson's Medical Supplies | |
| URBAN | Uptown Cleaners | 1610 Regent Street Sudbury, ON P3E 705 522 0400 | Robinson's Medical Supplies | |
| URBAN | Coniston Guardian | 24 Amanda St. Coniston, ON POM 1M0 705-694-5346 | Robinson's Medical Supplies | |
| URBAN | Wilson Pharmacy Pharmasave | 10 Serpentine St. Copper Cliff, ON POM 1N0 705-682-2841 | Robinson's Medical Supplies | |
| URBAN | Nickel Centre Pharmacy, Pharmasave | 3098 Falconbridge Hwy, Garson, ON P3L 1P5 705-693-2738 | Robinson's Medical Supplies | |
| URBAN | Val Est Pharmacy, Pharmasave | 3179 Hwy 69 N, Unit # 1 Val Caron, ON P3N 1S5 705-897-1867 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| URBAN | Walden Medicine Shoppe | 135 Regional Rd. 24 Walden Plaza 705-692-5868 | Robinson's Medical Supplies | |

Deliveries during extended hours: Only available for New Patients or Patient Resumes/ Hospital Discharges/End of Life Private-Pay Patient Deliveries: Robinson's to coordinate for Patients when requested

*Note: Vendor will contact client when supplies are ready for pickup at depot or to coordinate private-pay patient home delivery requests

SUDBURY-MANITOULIN-PARRY SOUND SUB-REGION

(Questions or issues call: Robinson's Pharmasave (705-566-5551 Option 4)

| Service Area | Depot/Drop-Off Location | Pick-Up Location | Service Type | Key Notes |
|-----------------|---|---|---|-----------|
| URBAN | Capreol Rx | 21 Young Street Capreol, ON P0M 1H0 705-858-8888 Fax: 705-858-8088 | Robinson's Medical Supplies | |
| URBAN | Dowling Pharmacy | 30 Main St. Dowling, ON 705-855-4562 | Robinson's Medical Supplies | |
| URBAN | Levack Pharmacy | 2 – 1 st Ave North Onaping Falls, ON 705-966-3612 | Robinson's Medical Supplies | |
| RURAL | Hagar General Store | 7147 Hwy 17 E, Hagar, ON 705-207-1006 | Robinson's Medical Supplies | |
| RURAL | Noelville Pharmacy | 27 St David N, Noelville, ON POO 2N0 705-898-2999 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| RURAL | St. Charles Cooperative Store | 25 Main St. St. Charles, ON 705-867-2000 | Robinson's Medical Supplies | |
| RURAL | Pitfield General Store | 7 Channel Killarney, ON P0M 2A0 705-287-2872 | Robinson's Medical Supplies | |
| RURAL | Robinsons Pharmasave | 119 Tudhope Espanola, ON 705-869-1561 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| RURAL | Home Hardware | 275 Imperial St Massey, ON P0P 1P0 705-865-2135 | Robinson's Medical Supplies | |
| RURAL | Little Current Guardian | 41 Meredith St Little Current, ON POP 1K0 705-368-2530 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| RURAL | Little Current Hospital Nursing Station | 11 Merideth Street E. Little Current, ON P0P 1K0 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| RURAL | Manitowaning Pharmacy | 1 Wellington Manitowaning, ON P0P 1N0 705-859-2323 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| RURAL | Mindemoya Pharmacy | Guardian Pharmacy 11 Yonge St. Mindemoya, ON P0P 1S0 705-377-5411 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| RURAL | Robertson IDA | 39 Meredith Street Gore Bay, ON POP 1H0 705-282-2147 | Robinson's Medical Supplies and IV Equipment & IV Supplies | |
| RURAL | Central Guardian Pharmacy | 45 Meredith Street Gore Bay, ON P0P 1H0 705-282-2409 | Robinson's Medical Supplies | |

Deliveries during extended hours: Only available for New Patients or Patient Resumes/ Hospital Discharges/End of Life **Private-Pay Patient Deliveries:** Robinson's to coordinate for Patients when requested

*Note: Vendor will contact client when supplies are ready for pickup at depot or to coordinate private-pay patient home delivery requests

SUDBURY-MANITOULIN-PARRY SOUND SUB-REGION

(Questions or issues call: Robinson's Pharmasave (705-566-5551 Option 4)

| Service Area | Depot/Drop-Off Location | Pick-Up Location | Service Type | Key Notes |
|-----------------|--|--|---|--|
| RURAL | Britt Nursing Station | Box 75 987 Riverside Dr. Britt, ON POG 1A0 705-383-2375 | Robinson's Medical Supplies | |
| RURAL | Duck Rock General Store | 2093 124 Hwy Dunchurch, ON P0A 1G0 705-389-3322 | Robinson's Medical Supplies | |
| RURAL | Lane Family Pharmasave | 4 Pine Parry Sound Parry Sound, ON P2A 3B8 705-746-2187 | Robinson's Medical Supplies & Infusion (IV) Equipment & Supplies | |
| RURAL | Rosseau General Store | 1 Rice Parry Sound, ON P2A 2W9 705-732-4479 | Robinson's Medical Supplies | |
| RURAL | Home and Community Care Support Services West Parry Sound Hospital | 6 Albert Street Parry Sound, ON P2A 3A4 | Robinson's Medical Supplies (not for IV) | For Home and Community Care Support Services staff only Not for Patient pick-up For Hospital Discharge only (to facilitate discharge) |

Deliveries during extended hours: Only available for New Patients or Patient Resumes/ Hospital Discharges/End of Life **Private-Pay Patient Deliveries:** Robinson's to coordinate for Patients when requested

HOME AND COMMUNITY CARE SUPPORT SERVICES North East

www.nelhin.on.ca

 SUDBURY OFFICE
 BUREAU DE SUDBURY

 40 Elm Street, Suite 41-C
 40, rue Elm, bureau 41-C

 Sudbury ON P3C 1S8
 Sudbury ON P3C 1S8

 Tel : 705 522-3461
 Téléphone : 705 522-3461

 Toll Free : 1 800 461-2919
 Sans frais : 1 800 461-2919

 Fax : 705 522-3855
 Télécopieur : 705 522-3855

 TTY: 711 (ask operator 1-888-533-2222)
 ATS : 711 (veuillez demander au téléphoniste le 1-888-533-2222)

